

MTL TR 89-36

AD

AD-A208 755

AUTOMATED IMMERSION TESTING ROBOTIC SYSTEM - HARDWARE AND SOFTWARE DESIGN

SUZANNE G. W. DUNN
POLYMER RESEARCH BRANCH

April 1989

Approved for public release; distribution unlimited.



US ARMY
LABORATORY COMMAND
MATERIALS TECHNOLOGY LABORATORY



U.S. ARMY MATERIALS TECHNOLOGY LABORATORY
Watertown, Massachusetts 02172-0001

SDTICD
ELECTE
JUN 08 1989
H

89 6 08 093

The findings in this report are not to be construed as an official Department of the Army position, unless so designated by other authorized documents.

Mention of any trade names or manufacturers in this report shall not be construed as advertising nor as an official indorsement or approval of such products or companies by the United States Government.

DISPOSITION INSTRUCTIONS

Destroy this report when it is no longer needed.
Do not return it to the originator.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER MTL TR 89-36	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) AUTOMATED IMMERSION TESTING ROBOTIC SYSTEM - HARDWARE AND SOFTWARE DESIGN		5. TYPE OF REPORT & PERIOD COVERED Final Report
7. AUTHOR(s) Suzanne G. W. Dunn		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS U.S. Army Materials Technology Laboratory Watertown, Massachusetts 02172-0001 SLCMT-EMP		8. CONTRACT OR GRANT NUMBER(s)
11. CONTROLLING OFFICE NAME AND ADDRESS U.S. Army Laboratory Command 2800 Powder Mill Road Adelphi, Maryland 20783-1145		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS D/A Project: 1L162105.AH84
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE April 1989
		13. NUMBER OF PAGES 59
		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
18. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number)		
Laboratory automation Immersion testing Robotics	Composites Materials durability Environmental testing	Polymers Liquid diffusion Moisture uptake
20. ABSTRACT (Continue on reverse side if necessary and identify by block number)		
(SEE REVERSE SIDE)		

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

Block No. 20

ABSTRACT

Hardware and software development of a fully automated robotic work cell to determine the water sorption-diffusion properties of plastics and organic matrix composite materials and to handle all data acquisition is described. Immersion testing involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. This test method is labor intensive and tedious, involving repetitive operations such as transferring liquids and specimens, blotting and weighing specimens, and recording data. Immersion studies can take years to complete and generate enormous amounts of data which must be recorded and evaluated. Through automation, immersion testing is no longer limited to an 8-hour work schedule, labor costs have been reduced, productivity has been increased tenfold, and the reliability, precision, and accuracy of the test have been improved.

CONTENTS

	Page
I. INTRODUCTION	1
II. ANALYTICAL PROCEDURE	1
Manual Test Method	2
Automated Test Method	2
III. ANALYSIS OF MOISTURE SORPTION DATA	2
IV. DESIGN OBJECTIVES AND CONSTRAINTS	4
System Constraints	4
Design Objectives	5
V. HARDWARE DESIGN	
Equipment	6
System Configuration	10
VI. SOFTWARE DESIGN	11
Method Controller	11
Manager	13
VII. SYSTEM VERIFICATION	
Testing and Verification	17
Productivity Enhancement	19
VIII. FUTURE ENHANCEMENTS	19
APPENDIX A. METHOD CONTROLLER PROGRAM CODE	22
APPENDIX B. MANAGER MODULE PROGRAM CODE	30



Accession For	
NTIS GPA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

I. INTRODUCTION

The moisture content of fiber-reinforced polymer matrix composite materials significantly affects their durability and mechanical properties. Immersion testing, used to investigate the moisture sorption-diffusion behavior of composites, involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. Immersion testing is tedious and time consuming. Depending upon the characteristics of the composite and temperature, a single test may run for several months, generating large amounts of data which must be collected and evaluated. Weight measurements are limited to an 8-hour work schedule, and are subject to operator errors. To successfully characterize the sorption-diffusion properties of composites, many immersion tests under various conditions are required. Manpower constraints limit the number and variety of conditions at which materials can be evaluated.

A fully automated robotic system has been developed for performing immersion studies in support of polymer research to increase the productivity for the evaluation of existing and next generation materials. Our goals in designing the automated immersion testing system were: (1) to emulate or improve upon manual operations, (2) to increase sample throughput, (3) to improve the precision and accuracy of measurements, and (4) to provide sufficient flexibility to handle different types of specimens and enable future expansion of operations.

State-of-the-art laboratory robotics coupled with today's personal computers provides a flexible means for addressing the problems associated with immersion testing procedures in the laboratory. In the following sections, this report describes the details of the immersion test method and the hardware and software design of the robotics work cell developed for performing automated immersion studies continuously and unattended. Section II (Analytical Procedure) describes the test method for performing manual and automated immersion studies. Section III (The Analysis of Moisture Sorption Data) discusses the method for interpreting immersion test data and how this information may be used in the evaluation of composite materials. Our objectives and design constraints, including the anticipated problems and concerns with maintaining the integrity of an immersion study using laboratory robotics, are discussed in Section IV (Design Objectives and Constraints). The hardware and software designs are detailed in Sections V (Hardware Design) and VI (Software Design), respectively. Section VII (System Verification) discusses the studies performed to compare the manual and automated test methods for validating the Automated Immersion Testing Robotic System, and illustrates the significant improvements achieved in the immersion test method with the use of laboratory robotics. Current plans for improving the design, and extending the immersion testing system to include accelerated aging studies in efforts to evaluate the durability and predict the lifecycle of composite materials, are discussed in Section VIII (Future Enhancements).

II. ANALYTICAL PROCEDURE

To fully understand the moisture sorption-diffusion behavior of one composite material processed under one set of processing parameters, numerous immersion tests under various temperature conditions are required. A single immersion test involves making periodic weight measurements and recording the change in the specimen's weight, with respect to the specimen's initial weight, as a function of immersion time. The following paragraphs describe the manual and automated immersion test method procedures, respectively.

Manual Test Method

Organic matrix composite test specimens are typically conditioned by drying under vacuum at 60°C (or below their glass transition temperature) for a period of 48 hours (or until there is no further decrease in weight). After allowing the specimens to equilibrate in a desiccator, their dimensions are measured and they are transferred into sample jars which contain 60 ml of the immersion liquid (often at elevated temperatures). Weight measurements are taken after the specimens have been immersed for 30 minutes, 1 hour, 2 hours, 4 hours, 6 hours, 24 hours, and 30 hours. Subsequent weight measurements are taken weekly until the specimens no longer exhibit a change in weight.

The weighing process involves removing the specimen from the immersion liquid, blotting the specimen with high absorbent filter paper to remove all surface liquid droplets, weighing the specimen, and returning it to the immersion liquid as quickly as possible. For each weight measurement, the total time the specimen has been immersed and the actual weight reading is recorded. By plotting the percent weight change versus the square root of time, the diffusion rate may be calculated using the Fickian law of diffusion, described in Section III (Analysis of Moisture Sorption Data).

Automated Test Method

In the initial designing phase of the immersion testing robotic work cell, the primary goal was to develop a system that performed immersion tests precisely the same way as in the established manual immersion testing procedures. Although the immersion test method is straightforward, many complex operations are required; e.g., manipulating small specimens, capping and uncapping specimen jars, blotting the surface liquid droplets from specimens, and making accurate weight measurements.

Initializing a specimen into the Automated Immersion Testing Robotic System involves manual as well as automated operations. First, specimens are conditioned manually (as stated in the Manual Test Method Section), their dimensional measurements are taken and entered into the computer by the operator, and the specimen is placed in a specimen holder within the robot's work envelope. The robot then proceeds by taking the specimen's "initial" weight measurement, placing the specimen in a designated specimen jar containing the immersion liquid at a given temperature, and returning the specimen jar to a temperature-controlled circulating water bath.

The process of making weight measurements for the remainder of the immersion test procedure follows the manual test method operations, as stated previously, with two major exceptions. First, the robot performs the required operations, and the computer records and calculates all immersion test data. Second, the robot is allowed to make weight measurements more frequently during the normal 8-hour work day, overnight, and on weekends.

III. ANALYSIS OF MOISTURE SORPTION DATA

The objective of the immersion test is to determine the moisture content (percent weight gain) of a material as a function of its immersion time. When water is used as the

immersion liquid, the test is essentially the ASTM Standard Test Method for Water Absorption of Plastics (D 570).¹ Determinations of the relative rate of water absorption are important in evaluating the effects of moisture exposure on such properties as mechanical strength, electrical resistivity, dielectric losses, physical dimensions, and appearance. Moisture content has significant effects on physical-mechanical properties and the long-term durability of polymeric and organic matrix composite materials.

To interpret immersion test data, moisture diffusion through the thickness of a test specimen can be described using a one-dimensional Fickian equation

$$dC/dt = D \, d^2C/dx^2 \quad (1)$$

where C is the moisture concentration in the specimen, x is the thickness parameter, t is exposure time, and D is the effective diffusivity of moisture through the thickness. Assuming that moisture diffusion and temperature are constant inside the specimen and under appropriate boundary conditions,² the solution for Equation 1 can be approximated as

$$C_i = G (C_e - C_o) + C_o \quad (2)$$

where C_i is the moisture content of the specimen at time t_i , C_e is the moisture content of the specimen at equilibrium (maximum moisture uptake), C_o is the initial moisture content (assumed to be uniform or zero), G is a time-dependent parameter defined as

$$G_i = 1 - \exp \left[-7.30 \left(\frac{Dt_i}{h^2} \right)^{0.75} \right], \quad (3)$$

and h is the thickness of the specimen exposed on both sides. Since (Equation 2) the percentage weight (wt%) gain $W_i = 100 G_i$, moisture diffusivity D_o can be determined directly from experimental data by plotting W_i versus $(t_i)^{1/2}$ and calculated either from the initial slope S_o

$$D_o = 0.196 \cdot h^2 \cdot (S_o/W_e)^2 \quad (4)$$

or from the time $t_{1/2}$ corresponding to one-half of the maximum moisture uptake value W_e for the specimen, as illustrated in Figure 1,

$$D_o = 0.049 \cdot (h^2/t_{1/2}) \quad (5)$$

The diffusivity D_o and equilibrium moisture content W_e are parameters which describe a specimen's moisture diffusion properties.

Figure 1 illustrates the typical sorption-diffusion curve following the Fickian law of diffusion.

1. ASTM D 570. *Standard Test Method for Water Absorption of Plastics*. 1987 Annual Book of ASTM Standards, v. 08.01, American Society for Testing and Materials, Philadelphia, PA, 1987.
2. SHEN, C.-W, and SPRINGER, G. S. *Moisture Absorption and Desorption of Composite Materials*. J. Composite Materials, v. 10, 1976, p. 2-20.




Figure 1. Typical moisture sorption/diffusion curve.

IV. DESIGN OBJECTIVES AND CONSTRAINTS

In addition to a flexible, programmable robot, the automation of immersion testing required peripheral equipment for: (1) maintaining specimen jars, immersion liquids, and test specimens at a constant and elevated temperature, (2) capping and uncapping specimen jars, (3) blotting specimens, and (4) making accurate weight measurements to the fifth decimal place. In the process of selecting and/or designing this specialized equipment, we found that the system could not be designed easily and reliably to handle all immersion tests in a generic sense. For this reason constraints were established.

System Constraints

The immersion testing robotic work cell was designed around two constraints; specimen configuration and limitations on immersion liquids.

Specimen Configuration: Since immersion specimens vary in size, shape, thickness, rigidity, and uniformity, standardized sample specifications had to be established. In this initial work, it was therefore decided to restrict the test to rigid, one-inch square or disk specimens of relatively uniform thickness (0.5 mm-1.9 mm). This decision was made for three reasons: (1) there is a high and continuing demand for evaluating composite materials; composite specimens are rigid and easily fabricated to the specified dimensions, (2) the specimen size is relatively consistent with ASTM test methods, and (3) uniform specimens of this configuration allow the robot the ability to handle specimens directly with custom forcepslike fingers, and allow the blotting operations to be simplified for enhanced reliability.

Immersion Liquids: Liquids used in immersion testing also vary. Some liquids are hazardous or corrosive, requiring specimens to be self-contained in closed jars. The majority of immersion tests, especially for testing composites, are performed with water or saline solutions. To provide as much flexibility as possible, a decision was made to demonstrate the

system using "safe" liquids, primarily distilled water, while simulating more hazardous situations by using screw-capped containers for future adaptations to more hazardous environments.

Design Objectives

As mentioned in the Introduction, our goals in designing the immersion testing system were: (1) to emulate or improve upon operations as specified in the manual immersion test method, (2) to increase sample throughput, (3) to improve the precision and accuracy of measurements, and (4) to provide sufficient flexibility to handle different types of specimens and enable future expansion of operations. Although there were many anticipated advantages to automating the immersion test procedure, there were also many concerns with maintaining the integrity of the test method. The following paragraphs address our design goals in greater depth, and point out advantageous as well as potential problem areas.

Improved Operations: Although the automated immersion test method duplicates the mechanics of the manual test method, several improvements were made in the test procedure. First, as discussed in Section II (Analytical Procedure), the manual test method limited weight measurements to five measurements on day 1, two measurements on day 2, one measurement on day 3, and one measurement per week for the remainder of the study. Through automation, weight measurements may be taken more frequently during a normal working schedule, and additional information may be obtained by making measurements overnight and during weekends and holidays. Using a simple algorithm, the robot can be programmed to look at past data to determine and alter the time interval between weight measurements as necessary.

One of the anticipated advantages to making more frequent weight measurements was to be able to capture information anomalous to Fickian's law of diffusion. Along with this expected advantage came a potential problem. Since the number of weight measurements taken during the course of an immersion study was significantly increased, we were concerned about the affects that this might have on the overall integrity of immersion test data (i.e., whether the increased amount of time that the specimens are removed from their respective liquid environments would affect overall test results).

The second major improvement in operations concerns specimen handling. In the manual test method, specimens are allowed to lay flat on the bottom of their respective sample jars, limiting the flow of immersion liquid around the specimens. To address this problem, and to provide a means for the robot to manipulate specimens into and out of specimen jars, sample jars were custom designed to hold rigid specimens vertically (refer to Section V, Hardware Design).

Sample Throughput: Although automation did not decrease the amount of time required to perform an immersion test, many more specimens may be evaluated during a period of time using robotics. Materials may now be evaluated under a variety of temperature conditions. Materials processed under varying processing conditions may now be studied. Along with performing the mechanical tasks of immersion testing, data handling, too, affects sample throughput. Robotics, coupled with computers, make the time required for data acquisition, data handling, and data reduction negligible.

Robotics generally include feedback and control systems which verify accuracy and functionality for each operation. For this reason, the automated system takes twice as long

to make a single weight measurement as a human operator. Determining what effects this additional length of time would have on the overall test result was a major concern.

Precision and Accuracy of Measurements: A large advantage of using robotics is that each weight measurement for each specimen is performed in exactly the same way. Measurements are taken exactly on time, and data acquisition is not subject to human error.

One major concern in maintaining and/or improving the accuracy of measurements was in the specimen blotting procedure. The manual blotting process relies on visual inspection to ensure that all water droplets have been removed from the specimen prior to a weight measurement (i.e., one droplet left on a specimen will result in a significant increase in weight). Since the automated system does not have that visual capability, the concern was in the reliability of removing all droplets from all surfaces, including specimen edges.

Flexibility: It was important to limit the scope of our design initially in order to produce a successful design with a limited cost investment. However, anticipating future applications, the robotic system was designed with the flexibility required to expand its capabilities and extend the variety of materials and test conditions.

V. HARDWARE DESIGN

Equipment

The Automated Immersion Testing Robotic System primarily consists of a host IBM PC/AT computer interfaced with a Zymark Zymate laboratory robotic work cell. In addition to a robot arm and controller, the robotic work cell is comprised of four major "work stations" designed for performing very specific tasks: (1) a temperature controlled, circulating water bath used to maintain the specimens and immersion liquids at a constant, and often elevated, temperature, (2) a jar capping/uncapping station, (3) a sample blotting station for removing surface droplets from the specimens prior to making weight measurements, and (4) a weighing station for measuring the weight of specimens at timed intervals. A description of each major work station and/or piece of equipment is given in the following paragraphs.

Robot and Controller: The robot and controller module include a Zymate I robot with a "general purpose" hand, a Zymate controller with an expansion memory card and EZLAB programming software, and a teaching pendant for remote positioning.

Robot Grippers: A Zymate general purpose hand was modified with two sets of custom fingers providing a dual capability. General grippers, made from aluminum with rubber padding on the inside, were designed to handle large sample jars from 1-3/4 to 3 inches in diameter. The specific dimensions are illustrated in Figure 2. Forcepslike fingers were designed to handle specimens directly (Figure 3). Also made of aluminum, these fingers are 5 inches in length and can handle a sample thickness of up to 1 inch.

Sample Jars: Standard off-the-shelf 100 mL (2 inches in diameter), Pyrex, screw-cap jars were individually modified to hold the immersion specimens in an upright, vertical position. The jars were modified in two ways: First, a custom, Pyrex specimen holder, consisting of two 1.5-inch concave pieces of glass, was fused 1.25 inches apart with a 2-mm slit opening to the inside bottom of each jar (Figure 4). The specimen holder is fluted at the top, allowing easy replacement of the specimen inside the jar. Second, a Pyrex "tab" was mounted inside

the 1.25-inch neck (inside diameter) of the jar, positioned perpendicular to the specimen holder. This tab protrudes out 0.125 inch towards the center of the specimen holder, and is used by the robot to locate the specimen. These modifications allow easy access to the specimens inside the jars, a more uniform flow of the immersion liquid around the specimens, and each specimen to have a unique liquid environment, if desired.




Figure 2. Robot gripper design.




Figure 3. Forceps finger design.




Figure 4. Specimen jar design.

Water Bath: A Fisher Scientific Model 260 temperature-controlled, circulating water bath is used to maintain the desired temperature of the specimens and their immersion liquids. The water bath was modified with a custom-designed rack consisting of an aluminum platform designed to fit within the inner chamber of the bath and a polypropylene platform with uniform jar slots which rests on top of the bath casing. Adjustable cylindrical fittings between the two platforms allow the sample jars to be immersed in the water bath at the appropriate depth.

The specimen jars are placed in the water bath so that only the tops of the jars are exposed. This allows the robot access to the specimen jars and a uniform temperature environment for exposure of specimens.

A blotter is used to remove excess water from the bottoms of sample jars as they are transferred from the water bath. This blotter consists of an absorbent sponge mounted on a polypropylene platform positioned 6 inches above the work bench.

Capping/Uncapping Station: A commercial capping station, Zymark Model Z410, is used for multiple tasks. Conventionally, it is used for capping and uncapping the sample jars. This "work station" is equipped with its own grippers, used to grip the sample jars securely, and a rotational base which rotates in the appropriate direction (clockwise for capping and counterclockwise for uncapping). Using the motion of the capping station, the robot, with a firm grip on the cap, lifts the cap up (or moves the cap downward, whichever the case may be) to complete the capping or uncapping procedures. A polypropylene cap holder, conforming to the shape of the inside of the jar caps, is mounted on a ring stand to provide a secure and accessible place for the robot to store the cap while attending to other tasks.

For our applications, the capping station is also used for "orienting" the jar, and thereby, locating the specimen in the jar. Using the capping station to rotate the jar, the robot's forcepslike fingers are partially inserted into the jar and opened. Contact between the fingers and the glass "tab," located in the neck of the jar, causes the jar to stop rotating and positions the robot fingers at a 90-degree angle directly above the specimen. In this way, the fingers can be lowered and always gain a firm purchase of the specimen for removal and handling.

Blotting Station: The blotting station (Figure 5) simulates the manual method of blotting specimens. The blotter consists of two 3- x 3.5-inch vertically mounted plates. The first plate, held by a ring stand, remains stationary, while the other plate, connected to an air-actuated pneumatic cylinder mounted on a ring stand, is allowed to move back and forth over a distance of 3 inches.




Figure 5. Specimen blotter design.

The design of the blotting station allows the robot's fingers to hold a specimen against a stationary plate while a mobile plate closes in on the specimen. Initially, the plates are positioned 2.75 inches apart. Highly absorbent filter paper is mounted on sponge pads located on the plates' inner surfaces. The filter paper collects surface droplets from the specimen and insures that the specimen edges are blotted as well. A relay switch, mounted to the mobile plate, is used to verify that the blotter is functioning properly at all times.

To insure removal of all surface droplets, each specimen must be regripped and blotted several times. An individual polypropylene specimen holder, mounted 6 inches above the working bench, allows the robot to set the specimen down and regrip it in a different location. This is done so that droplets that may have been trapped beneath the forcepslike fingers are removed. The filter paper, mounted on the mobile blotting plate, is extended 0.25 inch above the mobile plate to allow the robot to blot the inside of the forceps prior to regripping the specimen.

A heated air gun with a directional nose is positioned beneath the blotter to dry the filter paper after each specimen is blotted in preparation for the next specimen.

Weighing Station: A Mettler A163 dual-range, top-loading, analytical balance is used for weighing the specimens. A polypropylene specimen holder is mounted securely to the balance pan and positioned 1 inch below the top of the weighing cavity. Specimens are loaded into the top of the balance via the top sliding window. This window is automatically opened and closed using the Zymark Automatic Balance Door Opener mounted on top of the balance. The automatic door opener, which is connected to the window handle, slides the balance door back and forth using an air-actuated pneumatic cylinder. A relay switch is included for verification purposes.

The Zymark balance interface allows communication between the balance and the robot controller. This allows programmable control for taring the balance prior to each weight measurement.

Sensors: Two optical verification sensors (Model 37622) with LED and photocells are used to verify jar and cap manipulations. A third verification switch is used to verify that the specimens are retrieved from the various "work stations" properly.

Power and Event Controller: This module, Zymark Model Z830, allows programmable control of several logical inputs, switch closures, DC power supplies, and electrical power outlets. It is used for transmitting and receiving information to and from all of the various sensors, relay switches, and solenoid valves used by the Automated Immersion Testing Robotic System.

Host Computer: An IBM PC/AT with a 20 MB hard disk is used for data acquisition and handling. This computer includes a floppy diskette drive, DOS operating system, BASIC programming language, enhanced color monitor and adapter card, IBM Pro-printer, and RS1 graphics package.

The Zymark computer interface card, designed specifically for interfacing the Zymark controller with an IBM PC, allows real-time communication between the robot controller and the host computer.

System Configuration

The layout arrangement is illustrated in Figure 6. The robotic work cell is configured with each work station securely mounted on a 4 x 7 foot tabletop.




Figure 6. Immersion testing robotic system layout.

Hardware Interfacing/Connections: The robot controller is equipped with 16 communication ports which are necessary for interfacing the individual work stations and various modules into one complete system. Each standard Zymark "work station" and/or module includes an interfacing circuit board. A menu driven programming facility is provided with each interface board. Once these interfaces are connected to a communications port within the robot controller, the necessary programming tools for running and maintaining that particular module are installed. The robot arm, capping station, printer, and power and event controller modules only require a simple connection. The connection is made to the robot controller, via respective interfacing circuit boards and a communications port, to achieve software control over physical functions.

Host Computer: The host computer is interfaced with the robot controller, via a computer interface module connection, to one of the controller's communications ports. Communication between the two systems is achieved with a 9-pin to 25-pin RS232 connection between the controller's computer interface module card and the computer's RS232 communications port.

Pneumatic Cylinders: The blotting and weighing stations are powered by a 100 psi regulated air compressor connected to two 4-way solenoid valves (one for the blotter and one for the balance door opener). Vinyl tubing and a 3-way polypropylene connector connect the valves to the air compressor and the respective pneumatic cylinders. The auxiliary motion of

both work stations is programmable by a hardware connection between the solenoid valves and one of the eight switch closures, and the +12-volt power supply on the power and event controller (PEC) module. The status of the verification relay switch mounted on both stations is accessible through a connection between the relay switches and one of the PEC's eight logic inputs. Both of these relay switches are mounted so that the switch is "closed" when the pneumatic cylinders are compressed.

Balance: The communication between the balance (digital output) and the robot controller requires an RS232 interface board connection from the RS232 port of the balance to a communications port in the controller.

Air Dryer: Software control over the air dryer is achieved by plugging in the dryer to a 10-amp electrical outlet which is connected to one of the switch closures and the +5-volt power supply on the PEC.

Sensors: All three optical sensors are connected to a logical input on the PEC. This allows programmable access to the status of each sensor.

The placement of the two primary photo sensors is critical for detecting errors in the robot's routine. Since the majority of problems that may occur with the jar and cap manipulations center around the capping station, both sensors are focused on the capping station at a 90-degree angle to each other. Sensor 1 is positioned at a height required to detect several events: (1) if the cap is not capped properly (i.e., the cap is not threaded properly on the jar), (2) the cap is not removed from the jar during the uncapping procedure (i.e., the cap is fastened too tightly for the robot to release it from the jar), (3) the jar is lifted with the sample while trying to retrieve the specimen from the jar, (4) the jar is retrieved from the bath properly, and (5) the jar is retrieved from the capping station properly.

Sensor 2 is positioned at a height required to detect the following events: (1) the cap is removed from the jar during the uncapping procedure, (2) the cap is removed from the "cap holder" properly in order to recap the jar, (3) the jar is retrieved properly from the bath, and (4) the jar is retrieved properly from the capping station. The specifics of these error detection routines are discussed in greater detail in Section VI, Software Design.

VI. SOFTWARE DESIGN

The software design of the Automated Immersion Testing Robotic System is made up of two primary modules, the Method Controller and the Manager. The Method Controller is a set of subroutines written in the EZLAB programming language in the robot controller and used to instruct the robot through the mechanics of the immersion testing procedure. The Manager, written in BASIC on the IBM PC/AT, provides control of the sequence and scheduling of events, and handles all data acquisition. The Manager is the host program which calls on the Method Controller to instruct the robot through the appropriate set of tasks.

Method Controller

The main function of the Method Controller is to instruct the robot through the immersion testing procedure and handle the mechanics of the error detection, verification, and correction routines. All discrete robot movements are assigned and programmed using the

robot controller. Robot functions, or tasks, are programmed as subroutines using the EZLAB programming language. The organization of these tasks, which together make up the immersion test procedure, are coordinated in the Manager using BASIC on the IBM PC/AT.

A unique subroutine is created for each discrete task the robot performs. The primary subroutines which make up the immersion testing procedure of an immersion study are as follows:

1. remove sample jar from water bath,
2. blot bottom of jar (remove excess water from the bottom of the jar),
3. check for presence of jar (using two photo sensors, four positions on the jar are verified to ensure that the robot is handling the jar properly),
4. uncap the jar,
5. check for successful removal of cap,
6. place cap aside,
7. remove sample from jar,
8. check for successful removal of sample,
9. blot sample (remove surface droplets),
10. check for sample (after blotting procedure),
11. weigh sample (recording the weight),
12. check for sample removal (from balance),
13. place sample in jar,
14. cap sample jar,
15. check for presence of jar, and
16. return sample jar in bath.

Since the Manager makes the necessary decisions when procedural errors occur in the immersion test method routine, the above 16 primary subroutines are combined in higher level programs in the robot controller so that the last task in each program is a verification, or error checking, routine. When the Method Controller returns control to the Manager, the status of the sensor used for verification is also returned. This allows the Manager to make the appropriate decision and instruct the Method Controller to execute the most suitable action or correction subroutine. All subroutines and programs used to instruct the robot through the mechanics of the immersion test method are listed in Appendix A.

Manager

The Manager handles the scheduling activities, problem/error decisions, calculations, data storage and retrieval, and reporting features. The Manager, Appendix B, is divided into four primary sections: (1) the Initialization Procedure, (2) the Main Module or Scheduling Procedure, (3) the Immersion Testing Module, and (4) Function Key Subroutines. The Manager continuously searches for the next activity to be performed until one of the function keys, set up to terminate the program, is executed.

Initialization: The initialization portion of the program assigns all variables and flags, opens communication between the IBM and robot controller, opens all data files used to keep track of specimens, and initializes function keys (see section on function key subroutines). Six types of data files are used to keep track of specimens, events, and errors.

STUDIES.DAT is the data file which keeps track of all current specimens undergoing an immersion test. Each record in the data file contains information unique to a particular specimen: (1) the location (index) within the water bath; i.e., which jar contains the specimen, (2) the initial weight of the specimen, (3) the time the last weight measurement was made, (4) the date the last weight measurement was made, (5) the current time interval between weight measurements, (6) the time the next weight measurement should be made, and (7) the date the next weight measurement should be made.

JARSTATS.DAT is the data file which indicates the current status of each jar in the water bath (in this case, there are 30 bath locations indexed in a uniformed 6 x 5 rack). Each record (the record numbers correspond to the index of the jars within the bath) in this data file contains information unique to each specimen jar: (1) the rating of the jar (Since some jars are more reliable than others, they are rated so that specimens requiring more frequent weight measurements use the most reliable jars, and specimens weighed less often use less reliable jars. This limits the number of errors that occur while removing and replacing specimens in their immersion jars. A rating of "1" indicates the most reliable jars to be used for specimens which require measurements every 24 hours or less. A rating of "2" indicates the less reliable jars to be used for specimens which have a time interval between weighings of 25 or more hours.), (2) a flag to indicate if the jar is currently being used (1 indicates yes, 0 indicates no), (3) the current time interval between weighings assigned to the specimen in that particular jar, and (4) whether or not the jar contains the immersion liquid [1 indicates yes, 0 indicates no (empty)].

STATUS.REP is the data file which is for diagnostic purposes. As errors occur, a record of (1) the type of error, [0 indicates nonfatal, 1 indicates fatal (i.e., the study was discontinued for a particular specimen)], (2) which specimen/jar had the error, (3) the time of the error, and (4) a description of the error (what steps were taken to correct the error and if the error correction routines were successful). This data file is used by the operator in order to locate and fix, if appropriate, problems so that a study may be continued, and to locate problem jars (i.e., change their status) or problems in the routines which may require modifications.

DATAFILE.DAT is the data file which keeps track of file names used to store the raw data for each specimen. Each record of the data file corresponds to the unique specimen which occupies that record number (index) in the water bath. The name of the data files are assigned by the operator when a specimen is initialized.

The data files (unique for each specimen) contain the raw data for a unique specimen. Each record contains information related to one weight measurement: (1) the date of the measurement, (2) the time of the measurement, (3) the number of hours in the study past, (4) the weight measurement taken, and (5) the change in weight (current weight - initial weight).

SPECS.DAT is the data file which stores all information concerning the specimens currently undergoing the immersion test. Each record contains information unique to each specimen (the record number corresponds to the specimen/jar index number): (1) description of specimen, (2) type of resin, (3) type of fiber, (4) composite prepreg configuration (i.e., direction of fibers), (5) temperature of study, (6) specimen thickness, (7) specimen surface area, and (8) specimen conditioning parameters.

Main Module: The main module makes use of the STUDIES.DAT file by looping through all of the file's records. A comparison of the current date and time is made with the "next date" and "next time" fields of each record. If there is a match, or the current time and date is greater than the time and date of the next scheduled weighing, the immersion testing subroutine is invoked. After one pass through all records of the file (invoking the immersion testing subroutine as appropriate), this module checks to see if the operator has initiated any of the function keys (i.e., the operator, through the function keys, takes priority over the automation of the program).

There are seven function keys set up for the operator to: (1) stop the program and release control of the robot controller, (2) stop the program without releasing the controller, (3) plot/report the data on a particular specimen gathered so far, (4) start a new specimen, (5) initialize/change the status of the specimen jars in JARSTATS.DAT, or (6) print errors encountered thus far.

If no specimen at the current time needs a weight measurement, the robot is instructed to fill empty specimen jars with the immersion liquid (i.e., immersion liquids must reach the testing temperature before a specimen is initialized into the immersion study). The flow of this module is illustrated in the flow diagram in Figure 7.

The operator makes the decision to terminate the study on a particular specimen. In the event of termination, the operator terminates the Manager module, removes the specimen, uses function key 5 to change the status of that particular specimen jar, and restarts the Manager module. Once the jar's status has been altered appropriately, the system is designed to recognize the change and reuse the jar as needed.

Immersion Testing Subroutine: The immersion testing subroutine is separated into several categories or sets of lower level subroutines: (1) the primary subroutine for initiating a weight measurement on a specimen and all data acquisition, (2) the controller subroutines for communicating with the robot controller, setting variables, retrieving values of variables, and invoking EZLAB robot subroutines, (3) control routines for sequencing events for performing the immersion weight measurements (i.e., instructing the robot through the appropriate sequence of events and recognizing errors when they occur), and (4) error handling routines which are executed if errors are detected in the robot's procedure.

The abstract algorithm for making a weight measurement is as follows:

- The bath location (index) is set appropriately in the robot controller.
- The program to take a weight measurement is executed.
- The weight of the specimen is retrieved from the robot controller.
- The file name for data storage is retrieved from DATAFILE.DAT.
- The data file is opened and the current date, time, weight, calculated total number of hours in study, and calculated change in weight values are stored in the data file.
- The next weighing time and date are calculated and the time interval between weight measurements is increased, if appropriate. Currently, the algorithm for increasing the time interval takes the last five weight measurements of the study and averages them. If all five weight measurements are within ± 0.002 g of the mean, then the time interval between weight measurements is increased 100%. This information is then updated in STUDIES.DAT in the record associated with the specimen.
- Return to main program.




Figure 7. Immersion scheduling flow diagram.

Controller and Control Subroutines: The controller subroutines consist of repetitive portions of code primarily used for sending and retrieving information to and from the robot

controller. The control subroutines for running the immersion test are the actual system calls used to instruct the robot through the test procedure. This section of the code also includes all the error handling and decision sequencing.

Function Key Subroutines: Function keys 9 and 10 are used by the operator to terminate the Manager. If the operator needs to make changes in the Method Controller or reprogram robot positions, function key 10 must be used to release the control that the IBM has over the robot controller before the controller's keyboard will accept information. Function key 9 may be used when the operator does not require access to the robot controller.

Function keys 4, 5, and 6 perform their specific functions and return to the Manager's control. Function key 4 dumps the text found in STATUS.REP to the line printer and deletes the file. All information relating to errors is printed on the computer's line printer: (1) which jars or specimens encountered the errors, (2) what steps were taken to correct the problem, and (3) if the correction routines were successful. Function keys 5 and 6 are used to initialize and change, respectively, the flags which are used by the system to keep track of specimen jars.

The two most frequently used function key routines are function key 7 (starting a new specimen) and function key 8 (plotting data). The flow diagram, Figure 8, illustrates the algorithm for initializing a new specimen. Since some of the custom specimen jars are more reliable than others, jars are rated according to their reliability. In the process of initializing a new specimen to the system, each jar is evaluated sequentially in the following manner:

- If rating = 2 (less reliable), then next jar.
- If the jar is not in use and contains immersion liquid, then proceed with that jar.
- If the jar is not in use but does not contain the immersion liquid, then remember as a possible jar, but continue to next jar.
- If the jar is in use but time intervals between weighings is less than 24 hours, then remember as a possible jar for SWITCHING, but continue to next jar.

The most desirable situation is that a jar of a reliable rating is not being used and contains the immersion liquid. If such a jar is not available, then the next situation most desirable is a jar with a reliable rating which is empty, in which case the robot fills the jar, sets the timer for 1 hour to allow the immersion liquid to reach the desired temperature (i.e., that of the water bath), and moves on to other duties. When the hour is ended, the robot signals the operator and the initialization procedure continues.

If neither of the above two situations is true for any jar, then a jar which is currently used with a time interval greater than 24 hours is investigated. First, the system searches the status flags of each jar in JARSTATS.DAT looking for an unused jar with a rating equal to 2 (less reliable). If such a jar is found, the robot retrieves the first jar from the water bath (with a rating of 1 where time intervals are greater than 24 hours), removes the specimen, empties the jar, fills the jar with the immersion liquid, returns the jar to the water bath, and sets the timer for 1 hour. Next, the robot retrieves the second jar (with a rating of 2, and containing the immersion liquid), places the specimen in the jar, and returns the jar to the water bath. Again, the robot will move on to other tasks, as appropriate, until the timer indicates that 1 hour has passed. The robot signals the operator and the initialization procedure continues.




Figure 8. Initializing immersion study flow diagram.

Once the robot locates a jar for the new specimen, the system prompts the user for initial information and makes the appropriate changes in the files STUDIES.DAT, DATAFILE.DAT, JARSTATS.DAT, and SPECS.DAT.

The operator may use function key 8, at any time, to plot the data collected thus far for any specimen currently registered in the system. Function key 8 reformats the necessary information found in the specimen's unique data file, and uses the RS1 plotting language to plot the change in weight versus the square root of time (Figure 9 in Section VII, System Verification). Control is returned to the Manager at the operator's request.

Comments: The operator has the ability, through the use of the function keys, to override the Manager's control over the scheduling of events. Calculated times for the next weight measurements may not be precisely met; however, the times recorded for each weight measurement taken are the actual time of the measurements and not the calculated times.

VII. SYSTEM VERIFICATION

Testing and Verification

Two studies were performed to verify and test the reliability of the Automated Immersion Testing Robotic System. The first study involved performing accelerated immersion tests (the

immersion liquid was kept at a constant temperature of 60°C) on four glass-reinforced composite specimens machined from the same sheet of material. Two specimens were tested on the immersion system, while the other two specimens were tested by a human operator. This was done to compare the accuracy and to determine the robot's ability to duplicate, if not improve, the blotting and weighing routines. Weight measurements on the specimens were taken at the same time throughout a 3-month period, but not overnight or on weekends. The same water bath and sample jar types were used, and manual weight measurements were made with the same model of Mettler balance. Figures 9a and 9b show graphical results of specimens tested using the automated and manual methods, respectively. A comparison of individual data points indicates that the automated system can perform immersion studies with the same, if not better, accuracy as obtained by the manual method. Reliability was demonstrated in that the robotics system functioned autonomously during the test without committing an uncorrectable error. It is also apparent that the automated blotting procedure was performed reliably, since a single droplet of water left on the specimen during weighing produces a well-defined "outlier" data point and can be interpreted as such by the computer during analysis.




Figure 9. Immersion test data of an epoxy resin/aramid fiber laminate.

The second study involved performing accelerated immersion tests on four aramid-reinforced composite specimens. All four specimens were dried in a desiccator prior to the test. The only difference in this study was that the robot was allowed to run overnight and on weekends, while the operator kept to the same schedule as in the previous study. By increasing the number of weight measurements and making measurements during evenings and weekends, the Automated Immersion Testing Robotic System was able to demonstrate an

improvement in the precision and reliability of test results compared to the manual method as shown in Figure 9b.

Productivity Enhancement

A computer simulation was performed to illustrate the increase in productivity of the automated test method over the manual test method. The parameters used for each simulation are listed in Table 1. One of the advantages of the automated test method for performing immersion studies is that more weight measurements may be taken in the course of a single study. However, in order to compare the productivity of the two test methods, several constraints and assumptions were made. Since an immersion test takes the same number of days to complete a single sample independent of the test method, the automated method was constrained to the same number of weight measurements per specimen as the manual method. For this particular simulation, 15 weight measurements were taken in a period of 101 days for each specimen. It was assumed with the manual method that an operator works 8 hours per day and 47 weeks per year (i.e., 2 weeks of holidays, 2 weeks of vacation, and 1 week of sick leave) performing nothing but immersion testing tasks. For the manual method, it was also assumed that data handling (i.e., data entry and plotting/evaluating on a computer) is performed only after the 101-day study is completed. An operator requires 0.167 minute to enter the set of time and weight measurements in the computer, and 60 minutes to plot and evaluate each specimen. The automated method was allowed to perform 52 weeks per year, assuming no downtime, and the data handling time was assumed to be negligible. No constraints were put on the number of specimens that the automated system or manual system could handle at any time (i.e., it was assumed that there was an infinite number of specimen jars and bath slots). The manual blotting and weighing procedure takes approximately 2.5 minutes per specimen, while the automated blotting and weighing procedure takes approximately 6.0 minutes per specimen.

Table 1. SIMULATION PARAMETERS

Parameter	Manual Method	Automated Method
Time/Weight Measurement (min)	2.5	6.0
Working Hours/Day	8.0	24.0
Working Days/Week	5.0	7.0
Working Weeks/Year	47.0	52.0
Data Acquisition:		
Time/Specimen Plotted (min)	60.0	0.0
Data Entry (min)	0.167 * No. of Points	0.0

The simulation for each test method was performed to quantify the number of specimens that could be studied completely in one year. On the average, a human operator makes 29 weight measurements in an 8-hour day, while the Automated Immersion Testing Robotic System handles, on the average, 237 weight measurements in a 24-hour day. Thus, on a yearly basis, only 266 specimens can be tested manually while the automated system can test 2,775 specimens. Through automation, productivity has been increased by 1,170%.

VIII. FUTURE ENHANCEMENTS

Currently, plans for extending the Automated Immersion Testing Robotic System are underway. Figure 10 illustrates the design of an integrated system for materials durability

testing. This system will consist of four robot arms working simultaneously with two of the robots performing immersion testing as described in this report. Plans include the following additional capabilities: (1) automated specimen conditioning using a vacuum oven, (2) automated recycling of specimens, (3) a freeze-thaw study capability, (4) automated thickness measurements, (5) dry heat studies, (6) nondestructive testing, and (7) a variety of immersion temperature conditions. Enhanced software capabilities shall include: (1) sophisticated data analysis, (2) enhanced decision making capabilities through expert system development, and (3) machine learning capabilities.




Figure 10. Durability testing robotic system layout.

Robot A shall consist of: (1) a rack for holding unused specimen jars, (2) a station for immersion liquid disposal and liquid renewal, (3) a wash station for washing immersion specimens prior to weighings (e.g., if immersion liquid is a brine solution), (4) a specimen blotter, (5) a capping station, and (6) an analytical balance. Robot B shall consist of three fixed temperature block ovens set at 40°C, 60°C, and 80°C for maintaining immersion liquid temperatures and two capping stations.

To increase the number of weight measurements that can be performed in a 24-hour period, the immersion test procedure will be performed by dividing the test method tasks between robots A and B. Robot B shall be in charge of specimen jar manipulation into, and out of, the block ovens, and removing and replacing specimens in their respective containers. Robot A, once receiving a specimen from robot B, shall be in charge of blotting, washing, and weighing procedures. In addition to the immersion testing tasks, robot A shall also be in charge of jar preparations (i.e., cleaning, filling, and disposing of immersion liquids).

Specimens designated for accelerated aging studies shall be obtained from weathering chambers and environmental test sites, and will be placed on racks and introduced into the Durability Testing System at location C. Robot C shall be in charge of performing any further specimen-conditioning procedures, taking initial weight measurements, measuring specimen thicknesses, and handling specimen recycling/freezing-thaw procedures.

Robot D shall consist of variable temperature block ovens for performing dry-heat studies, and an NDT system for obtaining qualitative information on the moisture content and deterioration (delamination, crack formation, etc.) of immersion specimens. Several NDT methods are under investigation. Location D has been designed to allow system expansion. In the future, additional information may be obtained by incorporating dynamic mechanical analysis and chromatographic techniques for examining immersion specimens and liquids, respectively.

The Durability Testing System will allow numerous variations in specimen treatment/testing conditions. Specimens undergoing immersion testing may be dried and recycled, desorption measurements can be made, or the effects of freezing or high temperature treatment on NDT and sorption-diffusion properties may be studied. Complete records will be compiled for each specimen and a database will summarize results obtained for other materials/specimens. Depending upon relative changes in measurements and anticipated behavior, the robotics system may modify a specimen's test plan or alert the user of unusual behavior. When completed, the durability system will be integrated with a multitasking, multiprocessor, AI/Robotic work station.

REMOVE.JAR

```

get.jar.from.bath
blot.jar
check.for.jar

```

GET.JAR.FROM.BATH

```

gp.hand
move.to.bath.no.jar
over.cap.in.bath
grip.cap
raise.jar.from.bath
if bath.location<13 then 50
move.to.clearing
50

```

BLOT.JAR

```

move.to.blot.jar
lower.jar.to.sponge
set timer 1 3 seconds
wait for timer 1
move.to.blot.jar.2

```

CHECK.FOR.JAR

```

slant=0
jar.flag=0
sensor2.front
if photo.sensor.2=0 then 50
jar.flag=jar.flag+1
slant=1
50 sensor2.back
if photo.sensor.2=0 then 100
jar.flag=jar.flag+1
slant=2
100 move.to.sensor1
sensor1.front
if photo.sensor=0 then 150
jar.flag=jar.flag+1
slant=3
150 sensor1.back
if photo.sensor=0 then 200
jar.flag=jar.flag+1
slant=4
200 if slant=0 then 250
250

```

UNCAP.JAR

```

no module wait
move.to.capping.jar
module wait
open.capping.grippers
gp.hand
close.capping.grippers
repos.grippers.around.cat
grip.cap
uncap

```

UNCAP

```

no module wait
turn.to.uncap
module wait
robot.up.6mm
robot.up.8mm
robot.up.10mm
robot.up.12mm
lift.cap

```

CAP.ASIDE

```

avoid.sensor
place.cap.in.holding

```

PLACE.CAP.IN.HOLDING

```

move.to.cap.holding
lower.cap

```

gp.hand
move.to.cap.holding

GET.SAMPLE.FROM.JAR

move.to.capping.forceps
forceps.hand
orient.jar
ready.forceps.for.sample
close.capping.grippers
set timer 1 2 seconds
wait for timer 1
open.capping.grippers
lower.forceps.to.sample
grip.sample
set timer 1 2 seconds
wait for timer 1
lift.s
set timer 1 2 seconds
wait for timer 1
lift.sample.out.of.jar

ORIENT.JAR

lower.forceps
open.capping.grippers
open.forceps.orient
turn.to.orient

CHECK.FOR.SAMPLE.FROM.JAR

sample.flag=0
clear.beam.sensor
move.to.sensor.sample
set timer 1 2 seconds
wait for timer 1
if beam.sensor=0 then 30
sample.flag=1
30 clear.beam.sensor

GET.SPECIMEN

raise.arm.above.holding
forceps.hand
move.to.capping.forceps
orient.jar
forceps.hand
move.to.capping.forceps
dryer.on
blot.forceps.tips
move.towards.blotter
blot.forceps
dryer.off
to.sample.holder
get.sample.from.holder.1
check.for.sample.from.blotter

BLOT.FORCEPS.TIPS

forceps.hand
move.to.jar.blot.forceps
blot.tips
set timer 1 2 seconds
wait for timer 1
move.to.jar.blot.forceps

BLOT.FORCEPS

forceps.hand
blot.inside.forceps
move.in.blotter
set timer 1 10 seconds
wait for timer 1
move.towards.blotter

BLOT.INSIDE.FORCEPS

forceps.hand

GET.SAMPLE.FROM.HOLDER.1

sample.in.holder

```

move.to.sample.blotter
if blot.f.flag=1 then 5
move.to.left.blotter
move.down.to.blotter
goto 10
5 move.to.left.blotter.2
move.down.to.blotter.2
10 close.forceps
set timer 1 1 seconds
wait for timer 1
forceps.hand
if blot.f.flag=1 then 15
move.to.left.blotter
goto 20
15 move.to.left.blotter.2
20 move.to.sample.blotter

```

BLOT.SAMPLE

```

blot.f.flag=0
blot.flag=0
error.blot=0
blot.sample.center
if error.blot>0 then 100
put.sample.in.holder.1
blot.inside.forceps
get.sample.from.holder.2
blot.sample.left
if error.blot>0 then 100
put.sample.in.holder.2
blot.f.flag=1
blot.inside.forceps
blot.f.flag=0
get.sample.from.holder.3
blot.sample.right
if error.blot>0 then 100
put.sample.in.holder.3
get.sample.from.holder.1
100

```

OPEN.BLOTTER

```

error.blot=0
count=0
10 count=count+1
if count>5 then 50
blot.open
if blotter.switch=1 then 10
goto 100
50 error.blot=1
100

```

```

grip.sample
to.sample.holder.out

```

CHECK.FOR.SAMPLE.FROM.BLOTTER

```

sample.flag=0
move.to.sensor.sample
set timer 1 2 seconds
wait for timer 1
if beam.sensor=0 then 30
sample.flag=1
30 clear.beam.sensor

```

BLOT.SAMPLE.CENTER

```

hold.coount=1
move.to.sample.blotter
lower.sample.in.blotter
module wait
move.sample.right
blot.flag=1
close.blotter
if error.blot=1 then 50
set timer 1 1 seconds
wait for timer 1
blot.flag=2
open.blotter
if error.blot=1 then 50
lower.sample.in.blotter
50

```

CLOSE.BLOTTER

```

error.blot=0
count=0
5 count=count+1
if count>5 then 50
blot.close
if blotter.switch=0 then 5
goto 100
50 error.blot=1
100

```

PUT. SAMPLE. IN. HOLDER. 1

to.sample.holder
sample.in.holder
forceps.hand
to.sample.holder

GET. SAMPLE. FROM. HOLDER. 2

to.sample.holder.2
sample.in.holder.2
grip.sample
to.sample.holder.2

PUT. SAMPLE. IN. HOLDER. 2

to.sample.holder.2
sample.in.holder.2
forceps.hand
to.sample.holder.2

GET. SAMPLE. FROM. HOLDER. 3

to.sample.holder.3
sample.in.holder.3
grip.sample
to.sample.holder.3

PUT. SAMPLE. FROM. HOLDER. 3

to.sample.in.holder.3
sample.in.holder.3
forceps.hand
to.sample.holder.3

BLOT. SAMPLE. LEFT

hold.count=2
to.blotter.inside
inside.right
blot.flag=1
close.blotter
if error.blot=1 then 50
set timer 1 1 seconds
wait for timer 1
blot.flag=2
open.blotter
if error.blot=1 then 50
to.blotter.inside

50

BLOT. SAMPLE. RIGHT

hold.count=3
to.blotter.outside
outside.right
blot.flag=1
close.blotter
if error.blot=1 then 50
set timer 1 1 seconds
wait for timer 1
blot.flag=2
open.blotter
if error.blot=1 then 50
to.blotter.outside

50

TAKE. WEIGHT

tare.balance
weigh.sample

WEIGH. SAMPLE

if weight.type=1 then 1
move.to.balance.out

OPEN. BALANCE

count=0
5 count=count+1

```

    goto 2
1  move.to.balance
2  error.balance=0
   door.open.flag=1
   door.close.flag=0
   open.balance
   if error.balance>0 then 100
   if weight.type=1 then 3
   sample.in.balance.out
   goto 4
3  sample.in.balance
4  forceps.hand
7  move.to.balance.out
8  door.close.flag=1
   close.balance
   if error.balance>0 then 100
9  set timer 1 10 seconds
   wait for timer 1
   if weight.type=0 then 10
6  mw.sample=weight.sample
   goto 20
10 iw.sample(bath.location)=
   weight.sample
   if iw.sample(bath.location)=
   -999 then 20
20 door.open.flag=2
   open.balance
   if error.balance>0 then 100
   sample.in.balance.out
   grip.sample
   move.to.balance.out
   door.closed.flag=2
30 close.balance
100

```

CHECK. FOR. SAMPLE. FROM. BALANCE

```

sample.flag=0

clear.beam.sensor
move.to.sensor.sample
set timer 1 2 seconds
wait for timer 1
if beam.sensor=0 then 50
temp.weight=weight.sample
if temp.weight<0.1 then 45
pump.on
open.balance
move.to.balance
forceps.hand
sample.in.balance

grip.sample

```

```

if count>5 then 50
door.open
if door.status=1 then 5
goto 100
50 if door.open.flag=2 then 66
   error.balance=1
   goto 100
66 error.balance=2
100

```

CLOSE. BALANCE

```

count=0
5 count=count+1
if count>5 then 50
door.close
if door.status=0 then 5
goto 100
50 if door.closed.glag=2 then 75
   error.balance=3
   goto 100
75 error.balance=4
100

```

PUT. SAMPLE. IN. JAR

```

lift.sample.out.of.jar
lift.halfway
open.forceps.slightly
move.to.capping.forceps

```

GET. CAP

```

dryer.on
blot.forceps.tips
get.cap.from.holding
avoid.sensor
move.to.capping.cap

```

CAP. JAR

```

move.to.balance
pump.off
move.to.sensor.sample
set timer 1 2 seconds

```

```

cap.test=0
cap.sample.jar

```

```

wait for timer 1
if beam.sensor=0 then 50
45 sample.flag=1
stop.immersion.jar(bath.location)=1
50 clear.beam.sensor

```

GET.CAP.FROM.HOLDING

```

gp.hand
move.to.cap.holding
lower.cap.2
grip.cap
move.to.cap.holding

```

CAP.SAMPLE.JAR

GET.JAR

```

close.capping.grippers
move.to.capping.cap
place.cap.on.jar
no module wait
turn.capper.cap
module wait
robot.down.4mm
robot.down.6mm
robot.down.8mm
robot.down10mm
robot.down.12mm
set timer 1 3 seconds
wait for timer 1

```

```

open.capping.grippers
move.to.capping.jar
check.for.jar.2

```

CHECK.FOR.JAR.2

REPLACE.JAR

```

slant=0
jar.flag=0
sensor1.fron
if photo.sensor=0 then 50
jar.flag=jar.flag+1
slant=3
50 sensor1.back
if photo.sensor=0 then 100
jar.flag=jar.flag+1
slant=4
100 move.to.sensor1
sensor2.back
if photo.sensor.2=0 then 150
jar.flag=jar.flag+1
slant=2
150 sensor2.front
if photo.sensor.2=0 then 200
jar.flag=jar.flag+1
slant=1
200 move.to.blot.jar

```

```

move.to.blot.jar
dryer.off
put.sample.jar.in.bath
clear.all
gp.hand

```

PUT.SAMPLE.JAR.IN.BATH

```

if bath.location<13 then 10
move.to.clearing
10 move.to.bath.jar
lower.jar.in.bath
gp.hand
move.to.bath.no.jar

```


REGRIPPING

clear.sensor.2
regrip.jar
move.to.blot.jar
check.for.jar

PUT.JAR.IN.BATH

move.to.blot.jar
dryer.off
put.sample.jar.in.bath
clear.all
gp.hand

ERROR.BLOT.ROUTINE

lower.sample.in.blotter
move.to.sample.blotter

JAR.HANG.2

lift.halfway
grip.sample
close.capping.grippers
move.to.capping.forceps

GET.EXTRA.CAP

e.cap1
e.cap2
gp.wide
e.cap3
e.cap4
grip.cap
e.cap3
e.cap2
e.cap1

REGRIP.JAR

clear.all
box1
box2
gp.hand
box3
close.gp.hand
box4
box5
box6
box7
box8
box9
box10
box11
box12
box11
box4
box13
gp.hand
box14
grip.cap
box15
clear.all

NEXT.CAP

get.extra.cap
avoid.sensor
place.cap.in.holding
get.cap.from.holding
avoid.sensor
move.to.capping.cap

UNCAPPED.JAR.ASIDE

get.jar.from.capping.no.cap
clear.blot
put.jar.aside.no.cap

GET.JAR.FROM.CAPPING.NO.CAP

gp.hand
move.to.capping.to.grip.jar
move.to.grip.jar
open.capping.grippers
grip.jar
lift.jar.no.cap

PUT.JAR.ASIDE.NO.CAP

regrip1
regrip3
gp.hand
regrip1

CAP.ASIDE

avoid.sensor
place.cap.in.holding

FILL.JAR

w.pump.on
move.to.dispenser
valve.on
set timer 1 30 seconds
wait for timer 1
valve.off
w.pump.off
move.away.from.dispenser

MOVE.TO.DISPENSER

over.to.dispenser
in.to.dispenser
up.to.dispenser

MOVE.AWAY.FROM.DISPENSER

in.to.dispenser
over.to.dispenser

PUT.JAR.IN.CAPPING.NO.CAP

move.to.capping.to.grip.jar
open.capping.grippers
move.to.grip.jar
gp.hand
lift.cap
close.capping.grippers

SET.SAMPLE.ASIDE

clear.beam.sensor
blot.sample.center
put.sample.in.holder.1
gp.hand

EMPTY.JAR

move.to.beaker
mb.2
mb.3
mb.4
mb.5
mb.6
set timer 1 3 seconds
wait for timer 1
mb.4
mb.3
mb.2
move.to.beaker
gp.hand.jar

APPENDIX B. MANAGER MODULE PROGRAM CODE

```

80      ;Initialization of variables
100     OK$="OK"+CHR$(13)
120     NOTOK$="NOTOK"+CHR$(13)
140     ACK$="ACK"+CHR$(13)
160     LOOP=0
180     STARTFLAG=0
200     STARTBATH=0
220     STARTTIME$="00:00:00"
240     JAR.ASIDE.FLAG=0
280     DIM DAYS(12),TESTW$(5),JAR(30),LAST$(5)
300     DAYS(1)=31:DAYS(2)=28:DAYS(3)=31:DAYS(4)=30:DAYS(5)=31:DAYS(6)=30
320     DAYS(7)=31:DAYS(8)=31:DAYS(9)=30:DAYS(10)=31:DAYS(11)=30:DAYS(12)=31

340     ;Initialization of communication ports and data files
360     OPEN "COM1:1200,E,7,1,RS,DS0,LF" AS #1

380     OPEN "STUDIES.DAT" AS #2 LEN=53
400     FIELD #2, 3 AS BATH$, 7 AS IW$, 8 AS LTIME$, 10 AS LDATE$, 7 AS
        INTR$, 8 AS NTIME$, 10 AS NDATE$

        ;Data file STUDIES.DAT, length=53
        ;Purpose: Keeps track of current specimens undergoing immersion
        ;testing

        ;BATH$      Bath location          3      Integer
        ;IW$        Initial weight         7      Real
        ;LTIME$     Last weighing time     8      ###:##:##
        ;LDATE$     Last weighing date    10     ####:##:##
        ;INTR$      Current time interval  7      ###:##:##
        ;NTIME$     Next weighing time     8      ###:##:##
        ;NDATE$     Next weighing date    10     ####:##:##

420     OPEN "JARSTATS.DAT" AS #3 LEN=13
440     FIELD #3, 2 AS RATING$, 2 AS FU$, 7 AS CI$, 2 AS WWO$

        ;Data file JARSTATS.DAT, length=13
        ;Purpose: Status rating for all jars in water bath

        ;RATING$    1= 0 to 24 hr weighing intervals 2      Integer
        ;           2= 25 + hr weighing intervals
        ;FU$        jar in use/not in use             2      1/0
        ;CI$        current time interval             7      ####:##
        ;WWO$       jar with/with out water           2      1/0

460     OPEN "STATUS.REP" AS #5 LEN=128
480     FIELD #5, 2 AS TYPE$, 3 AS B$, 7 AS T$, 116 AS MESS$

        ;Data file STATUS.REP, length=128
        ;Purpose: Error diagnostics data file

        ;TYPE$      Type of error, fatal/non-fatal  2      0/1

```

;B\$	Bath location #	3	Integer
;T\$	Time of error	7	##:##:##
;MRSS\$	Error Message	116	String

```

500 LAST(2)=53
520 LAST(3)=13
540 LAST(5)=128
560 PRINT #1, "ENQ"+CHR$(13)
580 GOSUB 5680 ;Wait for controller to acknowledge

600 ;Initialization of Function keys
620 ON KEY(10) GOSUB 6480 ;Stop program and release robot control
640 ON KEY(9) GOSUB 6540 ;Stop program, do not release robot control
660 ON KEY(8) GOSUB 15380 ;Plot data
680 ON KEY(7) GOSUB 15760 ;Start new specimen
700 ON KEY(6) GOSUB 21640 ;Initialize data file JARSTATS.DAT
720 ON KEY(5) GOSUB 14520 ;Initialize/reinitialize Jar array
740 ON KEY(4) GOSUB 22940 ;Print status report
760 KEY(10) ON
780 KEY(9) ON
800 KEY(8) ON
820 KEY(7) ON
840 KEY(6) ON
860 KEY(5) ON
880 KEY(4) ON
900 KEY(10) OFF
920 KEY(9) OFF
940 KEY(8) OFF
960 KEY(7) OFF
980 KEY(6) OFF
1000 KEY(5) OFF
1020 KEY(4) OFF

1021 ;Function Key routines
1022 ;Key (10)—>Terminate program and release controller
1023 ;Key (9) —>Terminate program
1024 ;Key (8) —>Plotting/report generation
1025 ;Key (7) —>Start new specimen
1026 ;Key (6) —>Initialize JARSTATS.DAT
1027 ;Key (5) —>Initialize/Reinitialize Jar(?) array
1028 ;Key (4) —>Print Error Messages

```

;MAIN MODULE

```
1080 ;Determine if specimen(s) need weight measurements taken
1140 LOOP=LOOP+1
1160 ;Open STUDIES.DAT file and test next weight date and time of each
specimen
1220 IF LOF(2)=0 GOTO 1460 ;Test function keys
1240 FOR I=1 TO (LOF(2)/LAST(2))
1260     GET #2,I
1280     CDT$=DATE$ ;Current date
1300     CDATE$=(MID$(CDT$,7,4))+":"+ (MID$(CDT$,1,3))+":"+ (MID$(CDT$,4,2))
;Format Current date: Year:Month:Day
1320     CTM$=TIME$ ;Current time
1340     IF CDATE$>NDATE$ THEN GOSUB 1920 ;Take weight measurement
1360     IF (CDATE$=NDATE$) AND (CTM$>=NTIME$) THEN GOSUB 1920
1380 NEXT I

1420 ;Test function keys
1460 KEY(10) ON
1480 KEY(9) ON
1500 KEY(8) ON
1520 KEY(7) ON
1540 KEY(6) ON
1560 KEY(5) ON
1580 KEY(4) ON
1600 KEY(10) STOP
1620 KEY(9) STOP
1640 KEY(8) STOP
1660 KEY(7) STOP
1680 KEY(6) STOP
1700 KEY(5) STOP
1720 KEY(4) STOP

1740 ;Test for a new specimen waiting to be initialized
1760 ;Proceed with initialization if appropriate
1800 IF (STARTFLAG=1) AND (TIME$>STARTTIME$) THEN GOSUB 19640:LOOP=0:GOTO
1080
1820 IF LOOP=2 THEN GOSUB 14840:LOOP=0 ;If nothing to do, fill empty jars
1840 GOTO 1080 ;Restart Main Module
1860 CLOSE #1,#2,#3,#4,#5: END
```

;IMMERSION TEST SUBROUTINE

```

1920      ;Immersion Testing procedure
2040      ;Take weight measurement of specimen
2080      VARNAME$="BATH.LOCATION"
2100      BATH=VAL(BATH$)
2120      VARVALUE=BATH
2140      GOSUB 6000 ;Set BATH.LOCATION to BATH in controller
2180      VARNAME$="WEIGHT.TYPE"
2200      VARVALUE=1
2220      GOSUB 6000 ;Set WEIGHT.TYPE to 1 in controller
2260      IF JAR(BATH)=1 GOTO 1380 ;If error flag set, then move on to next jar
2280      GOSUB 6720 ;Take weight measurement
2300      VARNAME$="MW.SAMPLE"
2320      GOSUB 6360 ;Get value of MW.SAMPLE (weight of specimen)
2360      WEIGHT#=VARVALUE

2400      ;Locate specimen data file and update information
2440      OPEN "DATAFILE.DAT" AS #4 LEN=12
2460      FIELD #4, 12 AS FILE$

      ;Data file DATAFILE.DAT, length=12
      ;Purpose: Lists file names associated with each Jar in bath

      ;FILE$      File name and extension      12      String

2480      GET #4,BATH
2500      CLOSE #4

2540      ;Open data file
2580      OPEN FILE$ AS #4 LEN=48
2600      FIELD #4, 10 AS DDATE$, 8 AS DTIME$, 10 AS DHR$, 10 AS DWEIGHT$, 10
      AS DCHW$

      ;Data file *file*.DAT, length=48
      ;Purpose: Unique data file for each specimen

      ;DDATE$      Date of weighing      10      ####:###
      ;DTIME$      Time of weighing      8      ####:##
      ;DHR$        Total hours in study  10      real
      ;DWEIGHT$     Weight of specimen   10      real
      ;DCHW$        Change in weight     10      real

2620      GET #4,(LOF(4)/48)
2640      CHRS=VAL(DHR$) ;Hours in study
2660      CDT$=DATE$ ;Current date
2680      CDATE$=(MID$(CDT$,7,4))+ "-" +(MID$(CDT$,1,3))+ "-" +(MID$(CDT$,4,2))
      ;Reformat date
2700      CHW=WEIGHT#-VAL(IW$)

      ;Set date, time, weight and change in weight of specimen
2740      LSET DDATE$=CDT$
2760      LSET DTIME$=CTM$
2780      RSET DWEIGHT$=STR$(WEIGHT#)

```

```

2800      RSET DCHW$=STR$(CHW)

2820      ;Get last weighing time from file
2840      LH=VAL(MID$(LTIME$,1,2)) ;Hours
2860      LM=VAL(MID$(LTIME$,4,2)) ;Minutes
2880      LS=VAL(MID$(LTIME$,7,2)) ;Seconds

      ;Get current time
2900      CH=VAL(MID$(CTMS$,1,2))
2920      CM=VAL(MID$(CTMS$,4,2))
2940      CS=VAL(MID$(CTMS$,7,2))

2980      ;Calculate number of hours in current study
3020      IF CDATE$=LDATE$ GOTO 3800 ;If current day and the last day a
      weighing was made are equal

3040      ;Calculate number of hours in days past
3060      LD=VAL(MID$(LDATE$,9,2)) ;Last day
3080      CD=VAL(MID$(CDATE$,9,2)) ;Current day
3100      LMO=VAL(MID$(LDATE$,6,2)) ;Last month
3120      CMO=VAL(MID$(CDATE$,6,2)) ;Current month

3160      TH=23: TM=59: TS=59
3220      THRS=(TH-LH)+(TM-LM)/60+(TS-LS)/3600 ;Number of hours left in last
      day a weighing was made
3240      IF (MID$(CDATE$,1,4))=(MID$(LDATE$,1,4)) GOTO 3560 ;If the year has
      not changed
3280      THRS=THRS+24*(DAYS(LMO)-LD) ;Number of hours left in the last month a
      weighing was made
3300      TMO=LMO+1 ;next month

      ;Calculate number of hours left in the remainder of the year of the
      last weighing made
3320      WHILE TMO<=12
3340          THRS=THRS+24*(DAYS(TMO))
3360          TMO=TMO+1
3380      WEND

3400      ;Calculate number of hours in current year past
3420      TMO=1
3440      WHILE TMO<CMO ;While calculated month is less than current month
3460          THRS=THRS+24*(DAYS(TMO))
3480          TMO=TMO+1
3500      WEND

      ;Calculate number of hours past in current month
3520      THRS=THRS+24*(CD-1)
3540      GOTO 3940 ;Calculation finished

      ;Year has not changed
3560      IF CMO=LMO GOTO 3740
3580      THRS=THRS+24*(days(LMO)-LD) ;Hours left in month of last weighing
3600      TMO=LMO+1
3620      WHILE TMO<CMO

```

```

3640      THRS=THRS+24*(DAYS(TMO))
3660      TMO=TMO+1
3680      WEND
3700      THRS=THRS+24*(CD-1)
3720      GOTO 3940 ;Calculation finished

      ;Month has not changed
3740      THRS=THRS+24*(CD-LD-1) ;Days past in current month
3760      GOTO 3940 ;Calculation finished

3780      ;Day has not changed
3800      IF CS<LS THEN CS=CS+60: CM=CM-1
3820      IF CS>=LS THEN DS=CS-LS
3840      IF CM<LM THEN CM=CM+60: CH=CH-1
3860      IF CM>=LM THEN DM=CM-LM
3880      DH=CH-LH
3900      THRS=DH+DM/60+DS/3600 ;Calculation finished

3920      ;Store new information
3940      RSET DHR$=STR$(THRS+CHRS)
3960      PUT #4, (LOF(4)/48+1)

3980      ;Calculate next weighing time and date
4000      ;Increment time interval between weighings if appropriate
4060      IH=VAL(MID$(INSTR$,1,4))
4080      IM=VAL(MID$(INSTR$,6,2))
4100      NH=VAL(MID$(CTM$,1,2))
4120      NM=VAL(MID$(CTM$,4,2))
4140      NS=VAL(MID$(CTM$,7,2))
4160      NY=VAL(MID$(CDAT$,1,4))
4180      NMO=VAL(MID$(CDAT$,6,2))
4200      ND=VAL(MID$(CDAT$,9,2))
4220      CT=(LOF(4)/48-4) ;Looking at last five weight measurements
4240      IF CT<1 GOTO 4640 ;First weighing
4260      C=1
4280      TESTAV#=0
4300      FOR CC=CT TO (LOF(4)/48)
4320          GET #4, CC
4340          TESTW#(C)=VAL(DWEIGHT$)
4360          TESTAV#=TESTAV#+TESTW#(C)
4380          C=C+1
4400      NEXT CC
4420      AW#=TESTAV#/5
4440      FLAG=0
4460      FOR C=1 TO 5
4480          IF ABS(TESTAV#-TESTW#(C))>0.0002 THEN FLAG=FLAG+1
4500      NEXT C
4520      IF FLAG>1 GOTO 4640 ;Do not increase weight time increment
4540      TM=IM+IM:TH=0 ;Double time increment
4560      IF TM>=60 THEN TH=TH+1:TM=TM-60
4580      TH=TH+IH+IH
4600      IH=TH
4620      IM=TM
4640      GET #3,BATH

```



```

4660  TT$=STR$(IH):TT=LEN(TT$)-1:IH$=MID$(TT$,2,TT)
4680  ON TT GOSUB 4720,4740,4760,4770
4700  GOTO 4780
4720  IH$="000"+IH$: RETURN
4740  IH$="00"+IH$: RETURN
4760  IH$="0"+IH$: RETURN
4770  RETURN
4780  TT$=STR$(IM):TT=LEN(TT$)-1:IM$=MID$(TT$,2,TT)
4800  IF TT=1 THEN IM$="0"+IM$
4820  RSET CI$=IH$+" ":"+IM$
4840  PUT #3,BATH

4860  LSET LTIME$=CTM$
4880  LSET LDATE$=CDATE$
4900  RSET INTR$=IH$+" ":"+IM$
4920  TM=NM+IM: TH=0
4940  IF TM>=60 THEN TH=TH+1: TM=TM-60
4960  TH=TH+NH+IH
4980  IF TH<=23 THEN NH=TH: NM=TM: GOTO 5180

5000  WHILE (TH>23)AND(ND<DAYS(NMO))
5020      ND=ND+1
5040      TH=TH-24
5060  WEND
5080  IF TH<=23 GOTO 5180
5100  NMO=NMO+1
5120  ND=0
5140  IF NMO>12 THEN NY=NY+1:NMO=1
5160  GOTO 5000

5180  TT$=STR$(NY): TT=LEN(TT$)-1: NY$=MID$(TT$,2,TT)
5200  TT$=STR$(NMO): TT=LEN(TT$)-1: NMO$=MID$(TT$,2,TT)
5220  IF TT=1 THEN NMO$="0"+NMO$
5240  TT$=STR$(ND): TT=LEN(TT$)-1: ND$=MID$(TT$,2,TT)
5260  IF TT=1 THEN ND$="0"+ND$
5280  TT$=STR$(NH): TT=LEN(TT$)-1: NH$=MID$(TT$,2,TT)
5300  IF TT=1 THEN NH$="0"+NH$
5320  TT$=STR$(NM): TT=LEN(TT$)-1: NM$=MID$(TT$,2,TT)
5340  IF TT=1 THEN NM$="0"+NM$
5360  TT$=STR$(NS): TT=LEN(TT$)-1: NSS$=MID$(TT$,2,TT)
5380  IF TT=1 THEN NSS$="0"+NSS$
5400  LSET NDATE$+"-"+NMO$+"-"+ND$
5420  LSET NTIME$=NH$+" ":"+NM$+" ":"+NA$
5440  PUT #2,I
5460  CLOSE #4
5480  RETURN

```

;CONTROLLER SUBROUTINES

```
5640      ;Controller Ready?
5680      MESSAGE$=""
5700      WHILE (MESSAGE$<>ACK$)
5720          GOSUB 5840
5740      WEND
5760      RETURN

5800      ;Acquire messages from controller
5840      MESSAGE$=""
5860      CHARIN$=INPUT$(1,#1)
5880      IF CHARIN$=CHR$(10) THEN RETURN
5900      MESSAGE$=MESSAGE$+CHARIN$
5920      GOTO 5860

5960      ;Set the variable VARNAME to the value of VARVALUE in the controller
6000      PRINT #1, VARNAME$+"="+STR$(VARVALUE)+CHR$(13)
6020      GOSUB 6120
6040      RETURN

6080      ;Wait for response from controller
6120      MESSAGE$=""
6140      WHILE (MESSAGE$<>OK$)AND(MESSAGE$<>NOTOK$)
6160          GOSUB 5840
6180      WEND

6220      ;Test if communications problem
6260      IF MESSAGE$=CHR$(10)+"NOTOK" THEN RETURN 1860
6280      RETURN

6320      ;Retrieve value of VARNAME from controller
6360      PRINT #1, "? "+VARNAME$+CHR$(13)
6380      GOSUB 5840
6400      VARVALUE=VAL(MID$(MESSAGE$,3,LEN(MESSAGE$)-1))
6420      GOSUB 6120
6440      RETURN

6480      ;Stop program and disconnect communications
6500      PRINT #1, "BYE"+CHR$(13): CLOSE #1,#2,#3,#4,#5: RETURN 1860

6540      ;Stop program, do not disconnect communications
6560      CLOSE #1,#2,#3: RETURN 1860
```

;CONTROL SUBROUTINES FOR RUNNING DIMENSION TEST

```
6720      ;Take weight measurement
6760      GOSUB 7160 ;Get appropriate specimen jar, and uncap jar
6780      GOSUB 7480 ;Get specimen from jar
6800      GOSUB 7960 ;Blot specimen
6820      GOSUB 8260 ;Make weight measurement, return specimen to jar
6840      PROGNAMES$="DRYER.ON" ;dry blotter
6860      GOSUB 9060 ;Execute program
6880      GOSUB 8640 ;Cap jar and return to water bath
6900      RETURN

6940      ;Initialize a new specimen
6980      GOSUB 7140 ;Get appropriate specimen jar, and uncap jar
7000      GOSUB 7480 ;Get new specimen from specimen holder
7020      GOSUB 8260 ;Take initial weight measurement, return specimen to jar
7040      CTMS=TIMES ;Record initial time
7060      GOSUB 8640 ;Cap jar and return to water bath
7080      RETURN

7120      ;Get appropriate specimen jar, and uncap jar
7160      PROGNAMES$="REMOVE.JAR"
7180      GOSUB 9060 ;Execute program in controller
7200      VARNAME$="JAR.FLAG" ;Error flag for gripping jar properly
7220      GOSUB 6360 ;Get value of JAR.FLAG
7240      IF VARVALUE>0 THEN GOSUB 9260 ;Regrip jar properly
7260      PROGNAMES$="UNCAP.JAR"
7280      GOSUB 9060
7300      VARNAME$="PHOTO.SENSOR" ;Error flag for detecting jar
7320      GOSUB 6360
7340      IF VARVALUE=1 THEN GOSUB 10640 ;Cap not found, try uncapping again
7360      PROGNAMES$="CAP.ASIDE"
7380      GOSUB 9060
7400      RETURN

7440      ;Get specimen from jar
7480      PROGNAMES$="GET.SAMPLE.FROM.JAR"
7500      GOSUB 9060
7520      VARNAME$="PHOTO.SENSOR.2"
7540      GOSUB 6360
7560      IF VARVALUE=0 THEN GOSUB 11160 ;Jar lifted with specimen
7580      PROGNAMES$="CHECK.FOR.SAMPLE.FROM.JAR"
7600      GOSUB 9060
7620      VARNAME$="SAMPLE.FLAG"
7640      GOSUB 6360
7660      IF VARVALUE=1 THEN GOSUB 11380 ;sample not retrieved from jar
7680      CTMS=TIMES ;Record initial time of study
7700      RETURN

7740      ;Get new specimen from holder
7780      PROGNAMES$="GET.SPECIMEN"
7800      GOSUB 9060
7820      VARNAME$="SAMPLE.FLAG"
7840      GOSUB 6360
```

```

7860     IF VARVALUE=1 THEN GOSUB 14220 ;No specimen in holder
7880     RETURN

7920     ;Blot specimen
7960     CTM$=TIME$ ;Record time of weight measurement
7980     PROGNAMES$="BLOT.SAMPLE"
8000     GOSUB 9060
8020     VARNAME$="ERROR.BLOT"
8040     GOSUB 6360
8060     IF VARVALUE>0 THEN GOSUB 11840 ;Blotter malfunction
8080     PROGNAMES$="CHECK.FOR.SAMPLE.FROM.BLOTTER"
8100     GOSUB 9060
8120     VARNAME$="SAMPLE.FLAG"
8140     GOSUB 6360
8160     IF VARVALUE=1 THEN GOSUB 12640 ;Sample dropped during blotting
8180     RETURN

8220     ;Take weight measurement, return specimen to jar
8260     PROGNAMES$="TAKE.WEIGHT"
8280     GOSUB 9060
8300     VARNAME$="ERROR.BALANCE"
8320     GOSUB 6360
8340     IF VARVALUE>0 THEN GOSUB 12780 ;balance malfunction
8360     PROGNAMES$="CHECK.FOR.SAMPLE.FROM.BALANCE"
8380     GOSUB 9060
8400     VARNAME$="SAMPLE.FLAG"
8420     GOSUB 6360
8440     IF VARVALUE=1 THEN GOSUB 13360 ;Sample dropped in balance
8460     PROGNAMES$="PUT.SAMPLE.IN.JAR"
8480     GOSUB 9060
8500     VARNAME$="PHOTO.SENSOR.2"
8520     GOSUB 6360
8540     IF VARVALUE=0 THEN GOSUB 13520 ;jar lifted by forceps
8560     RETURN

8600     ;Cap jar and return to water bath
8640     PROGNAMES$="GET.CAP"
8660     GOSUB 9060
8680     VARNAME$="PHOTO.SENSOR"
8700     GOSUB 6360
8720     IF VARVALUE=1 THEN GOSUB 13660 ;cap not retrieved from holder
8740     PROGNAMES$="CAP.JAR"
8760     GOSUB 9060
8780     VARNAME$="PHOTO.SENSOR.2"
8800     GOSUB 6360
8820     IF VARVALUE=0 THEN GOSUB 13900 ;jar not capped properly
8840     PROGNAMES$="GET.JAR"
8860     GOSUB 9060
8880     VARNAME$="JAR.FLAG"
8900     GOSUB 6360
8920     IF VARVALUE>0 THEN GOSUB 9280 ;Regrip jar
8940     PROGNAMES$="REPLACE.JAR"
8960     GOSUB 9060
8980     RETURN

```

```
9020 ;Run program PROGNAMES$ in controller
9060 PRINT #1, PROGNAMES$+CHR$(13)
9080 GOSUB 6100
9100 RETURN
```

;ERROR HANDLING ROUTINES

```
9280 ;Jar not detected coming from water bath
9329 COUNT=0
9360 ;Try regripping jar no more than three times
9400 WHILE (VARVALUE>0) AND (COUNT<3)
9420     COUNT=COUNT+1
9440     MESSAGE$=STR$(VARVALUE)+" SIDE(S) OF SAMPLE JAR UNDETECTED."
9460     IF VARVALUE=4 GOTO 9920
9480     MESSAGE$=MESSAGE$+" REGRIPPING PROCEDURE EXICUTED.": GOSUB 10040
9520     RSET TYPE$="0": GOSUB 14360 ;Print message in STATUS.REP
9540 WEND
9580 IF VARVALUE=0 THEN RETURN ;Resume procedure if error corrected
9600 ON VARVALUE GOTO 9620,9800,9800,9920
9620 RSET TYPE$="1"
9640 MESSAGE$="REGRIPPING PROCEDURE UNSUCCESSFUL."

;reset error flag JAR.ASIDE.FLAG, flag specimen jar, print error
message, and move on to next jar
9680 IF JAR.ASIDE.FLAG=0 THEN MESSAGE$=MESSAGE$+" JAR SET ASIDE.":
JAR.ASIDE.FLAG=1: JAR(BATH)=1: GOSUB 14360: RETURN 1400

9720 MESSEAGE$=MESSAGE$+" NO ROOM TO PUT JAR ASIDE, SYSTEM SHUT DOWN
PROCEDURE EXICUTED."
9740 GOSUB 14360 ;Print message
9760 GOSUB 10200 ;System shut down procedure
9780 RETURN 1860 ;Terminate program

9800 RSET TYPE$=1
9820 MESSAGE$="JAR MAY NEED REPAINTING. REPLACING JAR IN BATH ATTEMPTED."
9840 GOSUB 10360 ;return jar to water bath
9860 JAR(BATH)=1 ;flag jar (i.e. discontinue study)
9880 GOSUB 14360 ;print message
9900 RETURN 1400 ;Next jar

9920 IF COUNT=1 THEN MESSAGE$=MESSAGE$+" JAR NOT PRESENT OR
UNRETRIEVABLE": RSET TYPE$="1": GOSUB 14360: JAR(BATH)=1: RETURN 1400

9960 MESSAGE$="JAR DROPPED DURING REGRIPPING PROCEDURE, SYSTEM SHUT DOWN
PROCEDURE EXICUTED.": GOSUB 14360: JAR(BATH)=1: GOSUB 10200

10000 RETURN 1860 ;Terminate program

10040 ;Regripping procedure
10080 PROGNAMES$="REGRIPPING"
10100 GOSUB 9060
10120 VARNAME$="JAR.FLAG"
10140 GOSUB 6360
10160 RETURN

10220 ;System shut down
10260 PROGNAMES$="DRYER.OFF"
10280 GOSUB 9060
```

```

10300 PRINT#1, "BYE"+CHR$(13)
10320 RETURN 1860

10360 ;Put jar in water bath
10400 PROGNAME$="PUT.JAR.IN.BATH"
10420 GOSUB 9060
10440 RETURN

10480 ;Get cap from holder
10520 PROGNAME$="AVOID.SENSOR": GOSUB 9060
10540 PROGNAME$="GET.CAP.FROM.HOLDING": GOSUB 9060
10560 PROGNAME$="AVOID.SENSOR": GOSUB 9060
10580 PROGNAME$="MOVE.TO.CAPPING.CAP": GOSUB 9060
10600 RETURN

10660 ;Cap not removed from jar properly
10700 COUNT=0
10720 WHILE (VARVALUE=1) AND (COUNT<3)
10740     COUNT=COUNT+1
10760     MESSAGE$="JAR UNCAPPED IMPROPERLY.  UNCAPPING PROCEDURE REPEATED."
10780     RSET TYPE$="0"
10800     GOSUB 14360
10820     PROGNAME$="GP.HAND": GOSUB 9060
10840     GOSUB 11000; Uncap jar
10860 WEND
10880 IF VARVALUE=0 THEN RETURN
10900 MESSAGE$="UNCAPPING UNSUCCESSFUL.  JAR REPLACED IN BATH."
10920 RSET TYPE$="1"
10940 GOSUB 14360; print message
10960 GOSUB 8840; Get jar, put jar in bath
10980 RETURN 1400 ;next jar

11020 ;Uncapping procedure
11060 PROGNAME$="UNCAP.JAR"
11080 GOSUB 9060
11100 VARNAME$="PHOTO.SENSOR"
11120 GOSUB 6360
11140 RETURN

11180 ;Jar lifted with sample
11220 MESSAGE$="JAR LIFTED WHILE TRYING TO RETRIEVE SPECIMEN FROM JAR.  JAR
REPLACED IN BATH.": GOSUB 14360
11240 RSET TYPE$="1"
11260 PROGNAME$="JAR.HAND1"
11280 GOSUB 9060
11300 GOSUB 8640; Get cap, cap jar, put jar in bath
11320 JAR(VAL(BATH$))=1
11340 RETURN 1400; Next jar

11380 ;sample not detected from jar
11420 MESSAGE$="SPECIMEN NOT RETECTED DURING REMOVAL ATTEMPT FROM JAR.
PROCEDURE REPEATED."
11440 RSET TYPE$="0"
11460 GOSUB 14360; Print message

```

```

11480      GOSUB 8460; Put sample in jar
11500      PROGNAMES$="GET.SAMPLE.FROM.JAR"
11520      GOSUB 9060
11540      PROGNAMES$="CHECK.FOR.SAMPLE.FROM.JAR"
11560      GOSUB 9060
11580      VARNAME$="SAMPLE.FLAG"
11600      GOSUB 6360
11620      IF VARVALUE=0 THEN RETURN
11640      MESSAGE$=MESSAGE$+"SPECIMEN RETRIEVAL UNSUCCESSFUL. JAR REPLACED IN
      BATH."
11660      RSET TYPE$="1":GOSUB 14360
11680      GOSUB 8460 ;Put sample in jar
11700      GOSUB 8640 ;Get cap, cap jar, put jar in bath
11720      JAR(BATH)=1
11740      RETURN 1400 ;Next jar

11800      ;Blotter malfunction
11840      MESSAGE$="BLOTTER MALFUNCTION: "
11860      ON VARVALUE GOSUB 12020,12060
11880      MESSAGE$=MESSAGE$+"SHUTDOWN PROCEDURE EXECUTED. "
11900      RSET TYPE$="1"
11920      GOSUB 14360 ;Print message
11940      ON VARVALUE GOTO 12180,1860
11960      RETURN

12020      MESSAGE$=MESSAGE$+"WILL NOT CLOSE. "
12040      RETURN

12060      MESSAGE$=MESSAGE$+"WILL NOT OPEN. "
12080      RETURN

12140      ;Determining type of malfunction
12180      VARNAME$="HOLD.COUNT"
12200      GOSUB 6360
12220      PROGNAMES$="ERROR.BLOT.ROUTINE"
12240      GOSUB 9060
12260      ON VARVALUE GOSUB 1244,12300,12360
12280      RETURN 1860

12300      PROGNAMES$="PUT.SAMPLE.IN.HOLDER.2"
12320      GOSUB 9060
12340      GOTO 12400

12360      PROGNAMES$="PUT.SAMPLE.IN.HOLDER.3"
12380      GOSUB 9060
12400      PROGNAMES$="GET.SAMPLE.FROM.HOLDER.1"
12420      GOSUB 9060
12440      GOSUB 8460 ;Put sample in jar
12460      GOSUB 8640; Get cap, cap jar, put jar in bath
12480      GOSUB 10200; System Shutdown
12500      RETURN

12560      ;Blotter will not close

```



```

12600    RETURN 1860; Terminate program

12660    ;Sample dropped during blotting
12700    MESSAGE$=MESSAGE$+"SAMPLE NOT DETECTED AFTER BLOTTING. JAR PUT BACK
        IN BATH. "
12720    RSET TYPE$="1": GOSUB 14360
12740    GOSUB 8640
12760    RETURN 1400

12800    ;Balance door malfunction
12840    MESSAGE$="BALANCE DOOR MALFUNCTION: "
12860    ON VARVALUE GOSUB 12960,13040,13140,13280
12880    RSET TYPE$="0"
12900    GOSUB 8640
12920    GOSUB 10220
12940    RETURN 1860

12980    MESSAGE$=MESSAGE$+"BALANCE DOOR WILL NOT OPEN; ERROR 1. SAMPLE AND
        JAR PUT BACK IN BATH"
13000    GOSUB 8460
13020    RETURN

13060    MESSAGE$=MESSAGE$+"BALANCE DOOR WILL NOT OPEN: ERROR 2. JAR PUT BACK
        IN BATH"
13080    PROGNAMES$="MOVE TO CAPPING.FORCEPS"
13100    GOSUB 9060
13120    RETURN

13160    MESSAGE$=MESSAGE$+"BALANCE DOOR WILL NOT CLOSE: ERROR 3. WEIGHT NOT
        RECORDED, SAMPLE AND JAR BUT BACK IN BATH"
13180    PROGNAMES$="SAMPLE.IN.BALANCE.OUT":GOSUB 9060
13200    PROGNAMES$="GRIP.SAMPLE":GOSUB 9060
13220    PROGNAMES$="MOVE.TO.BALANCE.OUT":GOSUB 9060
13240    GOSUB 8460
13260    RETURN

13300    MESSAGE$=MESSAGE$+"BALANCE DOOR WILL NOT CLOSE: ERROR 4. WEIGHT
        RECORDED, SAMPLE AND JAR RETURNED TO BATH"
13320    GOSUB 8460
13340    RETURN

13380    ;Sample dropped in balance
13420    MESSAGE$="SAMPLE DROPPED IN BALANCE."
13440    RSET TYPE$="1":GOSUB 14360
13460    PROGNAMES$="MOVE.TOWARDS.BLOTTER": GOSUB 20850
13480    GOSUB 8640
13500    RETURN 1400

13540    ;Jar lifted with forceps
13580    MESSAGE$="JAR LIFTED WITH FORCEPS AS SAMPLE IS PLACED IN JAR."
13600    RSET TYPE$="0":GOSUB 14360
13620    PROGNAMES$="JAR.HANG2":GOSUB 9060
13640    RETURN

```

```

13680 ;Cap missing from holder
13720 MESSAGE$="CAP MISSING FROM HOLDER. WILL GET EXTRA CAP IF AVAILABLE. "
13740 PROGNAME$="NEXT.CAP":GOSUB 9060
13760 VARNAME$="PHOTO.SENSOR":GOSUB 6360
13780 IF VARVALUE=0 THEN RSET TYPE$="0": GOSUB 14360:RETURN
13800 IF JAR.ASIDE.FLAG=1 THEN MESSAGE$=MESSAGE$+"NO EXTRA CAP FOUND. ASIDE
POSITION FULL. SYSTEM SHUT DOWN.":RSET TYPE$="1":GOSUB 14360:GOSUB
10260:RETURN 1860
13820 MESSAGE$=MESSAGE$+"NO EXTRA CAP FOUND, UNCAPPED JAR PUT ASIDE."
13840 RSET TYPE$="1":GOSUB 14360
13860 PROGNAME$="UNCAPPED.JAR.ASIDE":JAR.ASIDE.FLAG=1:GOSUB 9060
13880 RETURN 1400

13920 ;Jar not capped properly
13960 MESSAGE$="JAR NOT THREADED PROPERLY DURING CAPPING PROCEDURE.
UNCAPING/RECAPING ATTEMPTED."
13980 RSET TYPE$="0": GOSUB 14360
14000 GOSUB 11020
14020 IF VARVALUE=1 THEN GOSUB 10700; Jar not uncapped
14040 PROGNAME$="CAP.ASIDE":GOSUB 9060
14060 PROGNAME$="GET.CAP":GOSUB 9060
14080 VARNAME$="PHOTO.SENSOR":GOSUB 6360
14100 IF VARVALUE=0 THEN RETURN
14120 MESSAGE$="RECAPING PROCEDURE UNCSUCCESSFUL. "
14140 RSET TYPE$="1"
14160 IF JAR.ASIDE.FLAG=1 THEN MESSAGE$=MESSAGE$+"NO AVAILABLE POSITIONS TO
PLACE JAR ASIDE. SYSTEM SHUT DOWN.": GOSUB 14360:GOSUB 10220:RETURN
1860
14180 RETURN

14380 ;Print Message to STATUS.REP
14420 RSET B$=STR$(BATH)
14440 LSET T$=TIME$
14460 LSET MESS$+MESSAGE$
14480 PUT #5, (LOF(5)/128)+1
14500 RETURN

```

;FUNCTION KEY ROUTINES

```
14520 ;Initialize jar(?) Array (Function key 5)
14580 LOOP=0
14600 PRINT:PRINT "INITIALIZING STATUS OF SAMPLE JARS IN WATER BATH":PRINT
14620 INPUT "INITIALIE ALL JARS? (Y/N): ",ANS$
14640 IF ANS$="N" THEN GOTO 14740
14660 FOR C=1 TO 30 ;Number of jar locations in bath
14680     JAR(C)=0
14700 NEXT C
14720 RETURN

14740 INPUT "ENTER JAR #,VALUE # (0/1) $ENTER 0,0 TO END†: ",K,KV
14760 IF K=0 THEN RETURN
14780 JAR(K)=KV
14800 GOTO 14740

14840 ;Fill Empty jars
14880 FOR FILL=1 TO (LOF(3)/13)
14900     GET #3,FILL
14920     JR=VAL(RATINGS$)
14940     IF JR=0 GOTO 15020
14960     BATH=FILL
14980     JWOW=VAL(WWOS$)
15000     IF JWOW=0 THEN GOSUB 15060:RSET WWOS="1":PUT #3,BATH: GOTO 15040
15020 NEXT FILL
15040 RETURN

15060 ;Filling Jars
15120 VARNAME$="BATH.LOCATION"
15140 VARVALUE=BATH: GOSUB 6000
15160 GOSUB 7120; Get jar, uncap jar, put cap aside
15180 PROGNAME$="GET.JAR.FROM.CAPPING.NO.CAP":GOSUB 9060
15200 PROGNAME$="FILL.JAR":GOSUB 9060
15220 PROGNAME$="PUT.JAR.IN.CAPPING.NO.CAP":GOSUB 9060
15240 GOSUB 10460:GOSUB 8680; Cap jar, put jar in bath
15260 RETURN

15320 ;Plot data (Function key 8)
15380 INPUT "ENTER FILE TO BE PLOTTED ",FILES$
15400 CLOSE:OPEN FILES$ FOR INPUT AS #1
15420 OPEN "RS1.TMP" FOR OUTPUT AS #2
15440 FOR I=1 TO 6:INPUT #1, DUMMYZ: NEXT I
15460 INPUT #1, IWEIGHT
15480 FOR I=1 TO 6: INPUT #1,DUMMY$:NEXT I
15500 INPUT #1, DUM: INPUT #1,DUM
15520 WHILE NOT EOF(1)
15540     INPUT #1,ITIME,WEIGHT
15560     STIME=SQR(ITIME)
15580     SWEIGHT=100*((WEIGHT-IWEIGHT)/IWEIGHT)
15600     PRINT #2,USING "###.###"; SWEIGHT;
15620     PRINT #2,",";
15640     PRINT #2,STIME
15660 WEND
```

15680 LOOP=0
15700 CLOSE:RETURN 1860

;STARTING A NEW STUDY (OR SPECIMEN) (FUNCTION KEY 7)

```

15760 LOOP=0
15840 ONE=0:TWO=0
15880 ;Search for available jar with rating=1 (jar filled with water)
15900 FOR C=1 TO (LOF(3)/LAST(3))
15920   GET #3,C
15940   JR=VAL(RATING$)
15960   IF JR<>1 GOTO 16120 ;Don't consider jars with a rating=2
15980   JFU=VAL(FU$)
16000   JH=VAL(MID$(CI$,1,4))
16020   JM=VAL(MID$(CI$,6,2))
16040   JWOW=VAL(WWO$)
16060   IF (JFU=0) AND (JWOW=1) THEN BATH=C:GOTO 16260; Jar ready for use
16080   IF (JFU=0) AND (JWOW=0) THEN ONE=C:GOTO 16120; Possible jar,
      continue searching
16100   IF (JH>=24) AND (JAR(C)=0) THEN TWO=C; Possible jar for switching
16120 NEXT C
16140 IF ONE>0 THEN GOSUB 19080:GOTO 18940 ;Fill jar 'ONE' with water
16160 IF TWO>0 THEN GOSUB 19900:GOTO 18940 ;Switch jars
16180 PRINT:PRINT "NO JARS AVAILABLE AT THIS TIME...":RETURN

16220 ;User Input for specimen specifications
16260 GET #3,BATH
16280 INPUT "ENTER SAMPLE DESCRIPTION: ",DISC$
16300 INPUT "ENTER RESIN TYPE: ",RESIN$
16320 INPUT "ENTER FIBER TYPE AND NUMBER OF PLY: ",FIBER$
16340 INPUT "ENTER FIBER CONFIGURATION: ",CONF$
16360 INPUT "ENTER TEMPERATURE OF ENVIRONMENT: ",TEMP$
16380 INPUT "ENTER CONDITIONING PARAMETERS (25 CHAR): ",COND$
16400 INPUT "ENTER INITIAL TIME INTERVAL (MIN): ";INTERVAL
16420 PRINT:PRINT
16440 INPUT "ENTER FILE NAME FOR NEW SAMPLE: ",PFILE$
16460 PRINT:INPUT "ENTER SAMPLE THICKNESS: ",THICK$
16480 PRINT:INPUT "ENTER SAMPLE HEIGHT DIMENSION: ",A
16500 PRINT:INPUT "ENTER SAMPLE WIDTH DIMENSION: ",B
16520 PRINT:INPUT "ENTER SAMPLE DIAGONAL (RIGHT TOP TO LEFT BOTTOM): ",C
16540 PRINT:INPUT "ENTER SAMPLE DIAGONAL (LEFT TOP TO RIGHT BOTTOM): ",D
16560 JH=0:JM=INTERVAL
16580 IF INTERVAL >=60 THEN INTRZ=INTERVAL/60:JH=INTRZ:JM=INTERVAL-(JH*60)
16600 TT$=STR$(JH):TT=LEN(TT$)-1:JH$=MID$(TT$,2,TT)
16620 ON TT GOSUB 16660,16680,16700,16640
16640 GOTO 16720

16660 JH$="000"+JH$:RETURN
16680 JH$="00"+JH$:RETURN
16700 JH$="0"+JH$:RETURN

16720 TT$=STR$(JM):TT=LEN(TT$)-1:JM$=MID$(TT$,2,TT)
16740 IF TT=1 THEN JM$="0"+JM$
16760 RSET CI$=JH$+":"+JM$
16780 RSET FU$=STR$(1)
16800 RSET WWO$=STR$(1)
16820 PUT #3,BATH

```

```

16860 ;Calculate sample surface area
16920 AREA=A*B
16940 HYP=SQR(A^2+B^2)
16960 DIAG1=(HYP-C)/2
16980 DIAG2=(HYP-D)/2
17000 AREA1=4*DIAG1
17020 AREA2=4*DIAG2
17040 AREA=AREA-AREA1-AEA2

17120 ;Store initial data in files
17160 OPEN "DATAFILE.DAT" AS #4 LEN=12
17180 FIELD #4, 12 AS FILES
17200 LSET FILES=FFILES$
17220 PUT #4, BATH
17240 CLOSE #4

17280 RSET BATH$=STR$(BATH)
17300 RSET INTR$=JH$+"":JMS$
17320 PRINT:PRINT "PLACE SPECIMEN IN SAMPLE HOLDER"

17360 ;Take initial weighing
17400 PRINT:PRINT:PRINT "HIT ANY KEY TO CONTINUE"
17420 AS=INKEY$:IF AS="" THEN 17420
17440 VARNAME$="BATH.LOCATION"
17460 VARVALUE=BATH
17480 GOSUB 6000
17500 VARNAME$="WEIGHT.TYPE"
17520 VARVALUE=0
17540 GOSUB 6000
17560 GOSUB 6940;Take initial weight measurement

17580 ;Calculating and storing time for next interval
17620 CDT$=DATE$
17640 CDATE$=MID$(CDT$,7,4)+"-"+MID$(CDT$,1,3)+"-"+MID$(CDT$,4,2)
17660 LSET LTIME$=CTMS$
17680 LSET LDATE$=CDATE$
17700 VARNAME$="IW.SAMPLE(BATH.LOCATION)"
17720 GOSUB 6360
17740 IWEIGHT#=VARVALUE
17760 RSET IW#=STR$(IWEIGHT#)
17780 CH=VAL(MID$(CTMS$,1,2))
17800 CM=VAL(MID$(CTMS$,4,2))
17820 NH=0
17840 NM=CM+JM
17860 IF NM<60 THEN NM=NM-60:NH=NH+1
17880 NH=NH+CH+JH
17900 TT$=STR$(NH):TT=LEN(TT$)-1:NH$=MID$(TT$,2,TT)
17920 IF TT=1 THEN NH$="0"+NH$
17940 TT$=STR$(NM):TT=LEN(TT$)-1:NM$=MID$(TT$,2,TT)
17960 IF TTR=1 THEN NM$="0"+NM$
17980 IF NH<23 THEN LSET NTIME$=NH$+"":NM$+"":MID$(CTMS$,7,2):LSET
NDATE$=CDATE$:GOTO 18400
18000 NMO=VAL(MID$(CDT$,1,2))
18020 ND=VAL(MID$(CDT$,4,2))

```

```

18040  NY=VAL(MID$(CDT$,7,4))
18060  WHILE NH>23
18080      ND=ND+1
18100      IF ND>DAYS(NMO) THEN NMO=NMO+1:ND=1
18120      IF NMO>12 THEN NY=NY+1:NMO=1:ND=1
18140      NH=NH-24
18160  WEND
18180  TT$=STR$(NH):TT=LEN(TT$)-1:NH$=MID$(TT$,2,TT)
18200  IF TT=1 THEN NH$="0"+NH$
18220  TT$=STR$(NM):TT=LEN(TT$)-1:NM$=MID$(TT$,2,TT)
18240  IF TT=1 THEN NM$="0"+NM$
18260  TT$=STR$(NY):TT=LEN(TT$)-1:NY$=MID$(TT$,2,TT)
18280  TT$=STR$(NMO):TT=LEN(TT$)-1:NMO$=MID$(TT$,2,TT)
18300  IF TT=1 THEN NMO$="0"+NMO$
18320  TT$=STR$(ND):TT=LEN(TT$)-1:ND$=MID$(TT$,2,TT)
18340  IF TT=1 THEN ND$="0"+ND$
18360  LSET NTIME$=NH$+" ":"+NM$+" ":"+MID$(CTMS,7,2)
18380  LSET NDATE$=NY$+" ":"+NMO$+" ":"+ND$
18400  PUT #2, (LOF(2)/LAST(2))+1

18440  ;Creating datafile
18480  OPEN FFILE$ AS #4 LEN=48
18500  FILED #4, 10 AS DDATE$, 8 AS DTIME$, 10 AS DHRSS$, 10 AS DWEIGHT$, 10
    AS DCHW$
18520  LSET DDATE$=CDT$
18540  LSET DTIME$=CTMS
18560  RSET DHRSS$=STR$(0)
18580  RSET DWEIGHT$(STR$(IWEIGHT#)
18600  RSET DCHW$="0"
18620  PUT #4, LOF(4)/48+1
18640  CLOSE #4

18680  ;Storing sample specifications
18720  OPEN "SPECS.DAT: AS #4 LEN=125
18740  FIELD #4, 25 AS SD$, 15 AS SR$, 15 AS SF$, 15 AS SC$, 10 AS ST$, 10
    AS STH$, 10 AS SA$, 25 AS SCOND$
18760  LSET SD$=DISC$
18780  LSET SR$=RESIN$
18800  LSET SF$=FIBER$
18820  LSET SC$=CONF$
18840  LSET ST$=TEMP$
18860  LSET STH$=THICK$
18880  LSET SA$=STR$(AREA)
18900  LSET SCOND$=COND$
18920  PUT #4, BATH
18940  CLOSE #4
18960  RETURN

19000  ;Fill jar "ONE" with water, set startflag
19040  PRINT:PRINT "SPECIMEN JAR IS BEING PREPARED AT THE MOMENT, AT THE
    BEEP INTRODUCE THE SPECIMEN (1 HR.)"
19080  STARTBATH=ONE
19100  STARTFLAG=1
19120  VARNAME$="BATH.LOCATION"

```

```

19140  VARVALUE=STARTBATH
19160  GOSUB 6000
19180  GOSUB 7120; Get and uncap jar
19200  PROGNAMES$="GET.JAR.FROM.CAPPING.NO.CAP":GOSUB 9060
19220  PROGNAMES$="FILL.JAR":GOSUB 9060
19240  PROGNAMES$="PUT.JAR.IN.CAPPING.NO.CAP":GOSUB 9060
19260  GOSUB 10460:GOSUB 8680; Cap jar and return to bath
19280  CTMS=TIMES
19300  ST=VAL(MID$(CTM,1,2))
19320  ST=ST+1
19340  TT$=STR$(ST):TT=LEN(TT$)-1:STARTTIMES$=MID$(TT$,2,TT)
19360  IF TT=1 THEN STARTTIMES$="0"+STARTTIMES$
19380  STARTTIMES$=STARTTIMES$+MID$(CTMS,3,6)
19420  GET #3,STARTBATH
19440  RSET FUS$="1"
19460  RSET WWO$="1"
19480  PUT #3,STARTBATH
19500  RETURN

19660  ;Sounding BEEP for initializing specimen waiting
19680  SOUNDS$="GGGGGGGGGG"
19700  PLAY ">XSOUNDS$;"
19720  BATH=STARTBATH
19740  GOSUB 16260
19760  STARTFLAG=0
19780  STARTBATH=0
19800  STARTTIMES$="00:00:00"
19820  RETURN

19860  ;Search for jar with rating=2 for switching
19900  STARTBATH=TWO
19920  SWITCH=0
19960  ;Search for possible jar with rating=2 to switch sample with
19980  FOR C=1 TO (LOF(3)/LAST(3))
20000     GET #3,C
20020     JR=VAL(RATING$)
20040     IF JR<>2 GOTO 20120
20060     JFU=VAL(FUS$)
20080     JWO=VAL(WWO$)
20100     IF (JFU=0) AND (WWO=1) THEN SWITCH=C: GOTO 20240
20120  NEXT C
20140  PRINT "NO BATH LOCATIONS ARE AVAILABLE"
20160  STARTFLAG=0
20180  STARTBATH=0
20200  STARTTIMES$="00:00:00"
20220  RETURN

20260  ;Switch jars
20300  VARNAME$="BATH.LOCATION"
20320  VARVALUE=STARTBATH
20340  GOSUB 6000
20360  GOSUB 7120; Get and uncap jar
20380  GOSUB 7440; Get specimen from jar
20400  PROGNAMES$="SET.SAMPLE.ASIDE":GOSUB 9060

```



```

20420  PROGNAMES$="GET.JAR.FROM.CAPPING.NO.CAP":GOSUB 9060
20440  PROGNAMES$="EMPTY.JAR":GOSUB 9060
20460  PROGNAMES$="FILL.JAR":GOSUB 9060
20480  PROGNAMES$="PUT.JAR.IN.CAPPING.NO.CAP":GOSUB 9060
20500  GOSUB 10480:GOSUB 8680; Cap jar and return to bath
20520  STARTFLAG=1
20540  CTMS$=TIME$
20560  VARNAME$="BATH.LOCATION"
20580  VARVALUE=SWITCH
20600  GOSUB 6000
10620  GOSUB 7120; Get and uncap jar
20640  PROGNAMES$="GET.SAMPLE.FROM.JAR":GOSUB 9060; Orient jar
20660  GOSUB 7760; Get specimen from holder
20680  PROGNAMES$="CLEAR.BEAM.SENSOR":GOSUB 9060
20700  GOSUB 8460; Put sample in jar
20720  GOSUB 8600; Cap and return jar

20780  ;Update 'STUDIES.DAT'
20820  FOR J=1 TO LOF(2)/LAS(2)
20840      GET #2,J
20860      IF BATH=STARTBATH GOTO 20900
20880  NEXT J
20900  RSET BATH$=SWITCH
20920  PUT #2,J

20960  ;Update 'DATAFILE.DAT'
21000  OPEN "DATAFILE.DAT" AS #4 LEN=12
21020  FIELD#4, 12 AS FILE$
21040  GET #4, STARTBATH
21060  FFILE$=FILE$
21080  RSET FILE$="EMPTY"
21100  PUT #4,SWITCH
21120  RSET FILE$=FFILE$
21140  PUT #4,SWITCH
21160  CLOSE #4

21200  ;Update 'JARSTATS.DAT'
21240  GET #3,STARTBATH
21260  JCI$=CI$
21280  RSET CI$="0000:00"
21300  PUT #3,STARTBATH
21320  GET #3,SWITCH
21340  RSET FU$="1"
21360  RSET CI$=JCI$
21380  PUT #3,SWITCH

21420  ;Update SPECS.DAT
21460  OPEN "SPECS.DAT" AS #4 LEN=125
21480  FIELD #4, 25 AS SD$, 15 AS SR$, 15 AS SF$, 15 AS SC$, 10 AS ST$, 10
    AS STH$, 10 AS SA$, 25 AS SCND$
21500  GET #4,STARTBATH
21520  PUT #4,SWITCH
21540  GET #4,3
21560  PUT #4,STARTBATH

```

21580 CLOSE #4
21600 RETURN

```

21640 ;JARSTATS.DAT set-up (Function key 6)
21700 LOOP=0
21720 IF LOP(3)/13>0 THEN GOTO 22260
21740 FOR J=1 TO 11
21760     RSET RATINGS$="1"
21780     RSET FU$="0"
21800     LSET CI$="0000:00"
21820     RSET WWO$="0"
21840     PUT #3,J
21860 NEXT J
21880 GET #3,3
21900 RSET RATINGS$="0"
21920 PUT #3,3
21940 GET #3,4
21960 RSET RATINGS$="0"
21980 PUT #3,4
22000 FOR J=12 TO 20
22020     RSET RATINGS$="2"
22040     RSET FU$="0"
22060     RSET CI$="0000:00"
22080     RSET WWO$="0"
22100     PUT #3,J
22120 NEXT J
22120 FOR J=21 TO 30
22140     RSET RATINGS$="0"
22160     RSET FU$="0"
22180     RSET CI$="0000:00"
22200     RSET WWO$="0"
22240     PUT #3,J:NEXT J
22260 PRINT:PRINT "SET UP FILE JARSTATS.DAT"
22280 PRINT:PRINT "INFORMATION REQUIRED IS AS FOLLOS:"
22300 PRINT:PRINT "RATING: 0 - UNUSABLE"
22320 PRINT "      1 - < 24 HOURS"
22340 PRINT "      2 - >= 24 HOURS"
22360 PRINT:PRINT "USAGE: 1 - IN USE; 2 - NOT IN USE"
22380 PRINT:PRINT "TIME INTERVAL: XXXX:XX (HRS:MIN)"
22400 PRINT:PRINT "FILLED/NOT FILLED: 1-FILLED; 0-NOT FILLED"
22420 PRINT:INPUT :SETUP ALL BATH LOCATIONS? Y/N ",ANS$
22440 IF MID$(ANS$,1,1)="Y" GOTO 22780
22460 PRINT:INPUT "ENTER BATH LOCATION", ANS$
22480 IF ANS$="" GOTO 22920
22500 BATH=VAL(ANS$)
22520 GET #3,BATH
22540 PRINT:PRINT "RATING= ",RATINGS$
22560 PRINT:PRINT "USAGE= ",FU$
22580 PRINT:PRINT "INTERVAL= ",CI$
22600 PRINT:PRINT "FILLED/NOT FILLED= ",WWO$
22620 PRINT:INPUT "ENTER PARAMETERS (4): ",JR,JFU,JCI$,JWWO
22640 RSET RATINGS$=STR$(JR)
22660 RSET FU$=STR$(JFU)
22680 RSET CI$=JCI$
22700 RSET WWO$=STR$(JWWO)
22720 PUT #3,BATH
22740 GOTO 22460

```

```
22780   FOR J=1 TO 30
22800       PRINT:PRINT "ENTER WWO$ PARAMETER FOR BATH LOCATION ",J
22820       INPUT JR
22840       GET #3,J
22860       RSET WWO$=STR$(JR)
22880       PUT #3,J
22900   NEXT J
22920   RETURN
```

```

22940 ;Print status report and delete file (Function key 4)
23040 LPRINT " BATH TYPE OF ERROR"
23060 LPRINT " ERROR"
23080 LPRINT
23100 FOR ST=1 TO LOF(5)/LAST(5)
23120 GET #5,ST
23140 LPRINT " ",B$;
23160 LPRINT " ";
23180 TYPE=VAL(TYPE$)
23200 IF TYPE=0 THEN LPRINT "NOT FATAL ";
23220 IF TYPE=1 THEN LPRINT "FATAL ";
23240 VAL=LEN(MESS$)
23260 IF VAL<=50 THEN LPRINT MESS$:GOTO 23340
23280 MESSAGE$=MID$(MESS$,1,50):LPRINT MESSAGE$
23300 LPRINT " ";
23320 LPRINT MID$(,51,VAL-50)
23340 NEXT ST
23360 CLOSE #5
23380 OPEN "STATUS.REP" FOR OUTPUT AS #5
23400 CLOSE #5
23420 RETURN 460

```

DISTRIBUTION LIST

No. of Copies	To
1	Office of the Under Secretary of Defense for Research and Engineering, The Pentagon, Washington, DC 20301
	Commander, U.S. Army Laboratory Command, 2800 Powder Mill Road, Adelphi, MD 20783-1145
1	ATTN: AMSLC-IM-TL
	Commander, Defense Technical Information Center, Cameron Station, Building 5, 5010 Duke Street, Alexandria, VA 22304-6145
2	ATTN: DTIC-FDAC
1	Metals and Ceramics Information Center, Battelle Columbus Laboratories, 505 King Avenue, Columbus, OH 43201
	Commander, Army Research Office, P.O. Box 12211, Research Triangle Park, NC 27709-2211
1	ATTN: Information Processing Office
	Commander, U.S. Army Materiel Command, 5001 Eisenhower Avenue, Alexandria, VA 22333
1	ATTN: AMCLD
	Commander, U.S. Army Materiel Systems Analysis Activity, Aberdeen Proving Ground, MD 21005
1	ATTN: AMXSY-MP, H. Cohen
	Commander, U.S. Army Electronics Research and Development Command, Fort Monmouth, NJ 07703-5000
1	ATTN: AMDSD-L
1	AMDSD-E
	Commander, U.S. Army Missile Command, Redstone Scientific Information Center, Redstone Arsenal, AL 35898-5241
1	ATTN: AMSMI-RD-CS-R/Doc
1	AMSMI-RLM
	Commander, U.S. Army Armament, Munitions and Chemical Command, Dover, NJ 07801
2	ATTN: Technical Library
1	AMDAR-LCA, Mr. Harry E. Pebly, Jr., PLASTEC, Director
	Commander, U.S. Army Natick Research, Development, and Engineering Center, Natick, MA 01760
1	ATTN: Technical Library
	Commander, U.S. Army Satellite Communications Agency, Fort Monmouth, NJ 07703
1	ATTN: Technical Document Center
	Commander, U.S. Army Tank-Automotive Command, Warren, MI 48397-5000
1	ATTN: AMSTA-ZSK
2	AMSTA-TSL, Technical Library

No. of
Copies

To

1 Commander, White Sands Missile Range, NM 88002
1 ATTN: STEWS-WS-VT

1 President, Airborne, Electronics and Special Warfare Board, Fort Bragg,
NC 28307
1 ATTN: Library

1 Director, U.S. Army Ballistic Research Laboratory, Aberdeen Proving Ground,
MD 21005
1 ATTN: SLCBR-TSB-S (STINFO)

1 Commander, Dugway Proving Ground, Dugway, UT 84022
1 ATTN: Technical Library, Technical Information Division

1 Commander, Harry Diamond Laboratories, 2800 Powder Mill Road, Adelphi, MD 20783
1 ATTN: Technical Information Office

1 Director, Benet Weapons Laboratory, LCWSL, USA AMCCOM, Watervliet, NY 12189
1 ATTN: AMSMC-LCB-TL
1 AMSMC-LCB-R
1 AMSMC-LCB-RM
1 AMSMC-LCB-RP

1 Commander, U.S. Army Foreign Science and Technology Center, 220 7th Street, N.E.,
Charlottesville, VA 22901
1 ATTN: Military Tech

1 Commander, U.S. Army Aeromedical Research Unit, P.O. Box 577, Fort Rucker,
AL 36360
1 ATTN: Technical Library

1 Director, Eustis Directorate, U.S. Army Air Mobility Research and Development
Laboratory, Fort Eustis, VA 23604-5577
1 ATTN: SAVDL-E-MOS (AVSCOM)

1 U.S. Army Aviation Training Library, Fort Rucker, AL 36360
1 ATTN: Building 5906-5907

1 Commander, U.S. Army Agency for Aviation Safety, Fort Rucker, AL 36362
1 ATTN: Technical Library

1 Commander, USACDC Air Defense Agency, Fort Bliss, TX 79916
1 ATTN: Technical Library

1 Commander, U.S. Army Engineer School, Fort Belvoir, VA 22060
1 ATTN: Library

1 Commander, U.S. Army Engineer Waterways Experiment Station, P. O. Box 631,
Vicksburg, MS 39180
1 ATTN: Research Center Library

No. of
Copies

To

- Commandant, U.S. Army Quartermaster School, Fort Lee, VA 23801
1 ATTN: Quartermaster School Library
- Naval Research Laboratory, Washington, DC 20375
1 ATTN: Code 5830
2 Dr. G. R. Yoder - Code 6384
- Chief of Naval Research, Arlington, VA 22217
1 ATTN: Code 471
- 1 Edward J. Morrissey, AFWAL/MLTE, Wright-Patterson Air Force, Base, OH 45433
- Commander, U.S. Air Force Wright Aeronautical Laboratories,
Wright-Patterson Air Force Base, OH 45433
1 ATTN: AFWAL/MLC
1 AFWAL/MLLP, M. Forney, Jr.
1 AFWAL/MLBC, Mr. Stanley Schulman
- National Aeronautics and Space Administration, Marshall Space Flight Center,
Huntsville, AL 35812
1 ATTN: R. J. Schwinghammer, EH01, Dir, M&P Lab
1 Mr. W. A. Wilson, EH41, Bldg. 4612
- U.S. Department of Commerce, National Institute of Standards and Technology,
Gaithersburg, MD 20899
1 ATTN: Stephen M. Hsu, Chief, Ceramics Division, Institute for Materials
Science and Engineering
- 1 Committee on Marine Structures, Marine Board, National Research Council,
2101 Constitution Ave., N.W., Washington, DC 20418
- 1 Librarian, Materials Sciences Corporation, Guynedd Plaza 11, Bethlehem
Pike, Spring House, PA 19477
- 1 The Charles Stark Draper Laboratory, 68 Albany Street, Cambridge, MA 02139
- Wyman-Gordon Company, Worcester, MA 01601
1 ATTN: Technical Library
- Lockheed-Georgia Company, 86 South Cobb Drive, Marietta, GA 30063
1 ATTN: Materials and Processes Engineering Dept. 71-11, Zone 54
- General Dynamics, Convair Aerospace Division, P.O. Box 748, Fort Worth, TX 76101
1 ATTN: Mfg. Engineering Technical Library
- 1 Mechanical Properties Data Center, Belfour Stulen Inc., 13917 W. Bay Shore Drive,
Traverse City, MI 49684
- 1 Mr. R. J. Zentner, EAI Corporation, 626 Towne Center Drive, Suite 205,
Joppatowne, MD 21085-4440
- Director, U.S. Army Materials Technology Laboratory, Watertown, MA 02172-0001
2 ATTN: SLCMT-TML
1 Author

U.S. Army Materials Technology Laboratory
Watertown, Massachusetts 02172-0001
AUTOMATED IMMERSION TESTING ROBOTIC
SYSTEM - HARDWARE AND SOFTWARE DESIGN
Suzanne G. W. Dunn

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION

Key Words

Technical Report MTL TR 89-36, April 1989, 59 pp-
illus-table, D/A Project: 1L162105.AH84

Laboratory automation
Immersion testing
Robotics

Hardware and software development of a fully automated robotic work cell to determine the water sorption-diffusion properties of plastics and organic matrix composite materials and to handle all data acquisition is described. Immersion testing involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. This test method is labor intensive and tedious, involving repetitive operations such as transferring liquids and specimens, blotting and weighing specimens, and recording data. Immersion studies can take years to complete and generate enormous amounts of data which must be recorded and evaluated. Through automation, immersion testing is no longer limited to an 8-hour work schedule, labor costs have been reduced, productivity has been increased tenfold, and the reliability, precision, and accuracy of the test have been improved.

U.S. Army Materials Technology Laboratory
Watertown, Massachusetts 02172-0001
AUTOMATED IMMERSION TESTING ROBOTIC
SYSTEM - HARDWARE AND SOFTWARE DESIGN
Suzanne G. W. Dunn

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION

Key Words

Technical Report MTL TR 89-36, April 1989, 59 pp-
illus-table, D/A Project: 1L162105.AH84

Laboratory automation
Immersion testing
Robotics

Hardware and software development of a fully automated robotic work cell to determine the water sorption-diffusion properties of plastics and organic matrix composite materials and to handle all data acquisition is described. Immersion testing involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. This test method is labor intensive and tedious, involving repetitive operations such as transferring liquids and specimens, blotting and weighing specimens, and recording data. Immersion studies can take years to complete and generate enormous amounts of data which must be recorded and evaluated. Through automation, immersion testing is no longer limited to an 8-hour work schedule, labor costs have been reduced, productivity has been increased tenfold, and the reliability, precision, and accuracy of the test have been improved.

U.S. Army Materials Technology Laboratory
Watertown, Massachusetts 02172-0001
AUTOMATED IMMERSION TESTING ROBOTIC
SYSTEM - HARDWARE AND SOFTWARE DESIGN
Suzanne G. W. Dunn

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION

Key Words

Technical Report MTL TR 89-36, April 1989, 59 pp-
illus-table, D/A Project: 1L162105.AH84

Laboratory automation
Immersion testing
Robotics

Hardware and software development of a fully automated robotic work cell to determine the water sorption-diffusion properties of plastics and organic matrix composite materials and to handle all data acquisition is described. Immersion testing involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. This test method is labor intensive and tedious, involving repetitive operations such as transferring liquids and specimens, blotting and weighing specimens, and recording data. Immersion studies can take years to complete and generate enormous amounts of data which must be recorded and evaluated. Through automation, immersion testing is no longer limited to an 8-hour work schedule, labor costs have been reduced, productivity has been increased tenfold, and the reliability, precision, and accuracy of the test have been improved.

U.S. Army Materials Technology Laboratory
Watertown, Massachusetts 02172-0001
AUTOMATED IMMERSION TESTING ROBOTIC
SYSTEM - HARDWARE AND SOFTWARE DESIGN
Suzanne G. W. Dunn

AD UNCLASSIFIED
UNLIMITED DISTRIBUTION

Key Words

Technical Report MTL TR 89-36, April 1989, 59 pp-
illus-table, D/A Project: 1L162105.AH84

Laboratory automation
Immersion testing
Robotics

Hardware and software development of a fully automated robotic work cell to determine the water sorption-diffusion properties of plastics and organic matrix composite materials and to handle all data acquisition is described. Immersion testing involves measuring weight changes of specimens as a function of the time they are immersed in water at a particular temperature. This test method is labor intensive and tedious, involving repetitive operations such as transferring liquids and specimens, blotting and weighing specimens, and recording data. Immersion studies can take years to complete and generate enormous amounts of data which must be recorded and evaluated. Through automation, immersion testing is no longer limited to an 8-hour work schedule, labor costs have been reduced, productivity has been increased tenfold, and the reliability, precision, and accuracy of the test have been improved.